

Invacare® | *Softform Premier Active 2*

NEW



THREE TIER CARE

Warranty

- Foam 8 years
- Cover 4 years
- Pump 2 years
- Insert 2 years

New and improved static dynamic pressure reducing mattress

The clinically proven **Softform Premier Active** has undergone some changes and now incorporates a range of new features to enhance product performance and patient comfort.

The mattress retains its properties as a 'High/ Very High Risk' static surface, but should a patient require stepping up to a dynamic surface, a discreet pump can be fitted to the air insert, transforming the mattress into an alternating surface, delivering additional levels of pressure relief. There is no requirement to transfer the client, thereby minimising any discomfort and reducing manual handling. In addition, there is no need to wait for an alternative mattress to be ordered and delivered, meaning clinical intervention can start immediately.

The mattress is appropriate for patients in Acute, Community and Residential environments.



Yes, you can.®



Static Dynamic Mattress

The air insert, contained within the mattress allows stepping up and down between static and dynamic alternating systems with minimum patient disruption.

The cells alternate on a ten minute cycle beneath the foam insert, delivering a very low impact at the point of change. This patented technology provides the comfort of a static foam mattress with the therapeutic benefit of alternating pressure.

Multi-stretch vapour-permeable cover

The mattress cover is multi-stretch and conforms to the deep contoured castellations on the upper surface of

the mattress, minimising "rucking", and further reducing the potential for skin tissue damage.

How to 'step up' the mattress



Detach the pipes from the foot end of the mattress.



Fit pipes to the Active pump.



Simply switch on, no additional adjustments required.



The patient is now lying on an alternating surface.



Features and Options



New air cell design

The cells offer a stable support surface along with improved comfort for the patient. Advanced material used for the insert allows the foam to glide over the cells with ease, reducing friction.



Individual replaceable air cells

In case of accidental puncturing, individual air inserts can easily be replaced at ward level.



Channels in Blue foam U core

The tubes that supply air to the cells sit into securely designed channels, providing added comfort for the patient.




Cost saving modular construction

The top cover, bottom cover, foam insert, Blue U core and air cells are all easily replaceable.



Concealed pump hose

The pipes are fully covered with a welded sleeve to promote effective infection control.



High Risk static mattress

In the event of a power failure, the patient is still resting on a 'High Risk' pressure reducing surface, minimising the risk of pressure damage.



Cable and Hose management

If pump is in use, cables can be lifted off the floor and attached to the mattress. If not in use, the air hose can be fastened out of the way.




Technical data




Softform Premier Active 2

Width:
830 - 900 mm
(10 mm increments)



Length:
1810 - 2100 mm
(10 mm increments)



Height:
152 mm



Weight:
14 kg



Max. User Weight:
247.6 kg (39 stone)



Cleaning: 80°

Pump

Noise Level: 32Db
Cycle Time: 10 minutes
(1 in 2 cycle)

Mattress is manufactured to:

- European Standard: BS EN 597-1 Cig
- European Standard: BS EN 597-2 Match
- BS7177:1996 Crib 5.

Colour



Baltic Blue

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be
correct at time of print.
Invacare reserves the right to alter
product specifications without prior
consultation.

Softform Premier Active 2 - UK -
11/2011 - 1552287