# Invacare® Matrix®

Vi / Kidabra Vi

en	<b>Cushion</b>	User Manual . . . . .	1
de	<b>Kissen</b>	Gebrauchsanweisung . . . . .	3
fr	<b>Coussin</b>	Manuel d'utilisation . . . . .	6
it	<b>Cuscino</b>	Manuale d'uso . . . . .	9
es	<b>Cojín</b>	Manual del usuario . . . . .	11
pt	<b>Almofada</b>	Manual de utilização . . . . .	14
nl	<b>Kussen</b>	Gebruiksaanwijzing . . . . .	16
da	<b>Pude</b>	Brugsanvisning . . . . .	19
no	<b>Setepute</b>	Bruksanvisning . . . . .	21
sv	<b>Dyna</b>	Bruksanvisning . . . . .	24
fi	<b>Pehmuste</b>	Käyttöohje . . . . .	26

en

## I General

### I.1 General information

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions. For further information contact Invacare® in your country (addresses see back page of this manual).

### I.2 Symbols in this user manual

In this User Manual warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

 **WARNING**  
Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

 **CAUTION**  
Indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

 **IMPORTANT**  
Indicates a hazardous situation that could result in damage to property if it is not avoided.

 Gives useful tips, recommendations and information for efficient, trouble-free use.

 This product complies with the directive 93/42/EEC for medical products. The launch date for this product is specified in the CE declaration of conformity.

 Manufacturer

 This manual **MUST** be given to the user of the product. **BEFORE** using this product, read this manual and save for future reference.

### I.3 Warranty

We provide a manufacturer's warranty of two years for the product in accordance with our General Terms and Conditions of Business in the respective countries. Guarantee claims can only be made through the provider from whom the appliance was obtained.

#### Quality and flame retardancy

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485 and ISO 9001.

The Invacare® Matrix® Vi Cushion features the CE mark, in compliance with the Medical Device Directive 93/42/EEC Class I. Invacare® is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum. We use only REACH compliant materials and components.

The foam and cover used to manufacture the Invacare® Matrix® Vi Cushion are fire safety tested and certified in accordance with EN 1021-1 & -2.

For further information please contact Invacare® in your country (addresses see back page of this manual).

### I.4 Intended use

This product has been designed to deliver effective pressure reduction and postural support to users. Normal use is defined by Invacare Ltd as when the product is placed without additional covers or padding between the user and the support surface. Replace the cover in case of damage resulting in reduced tightness or elasticity of the protective barrier.

### I.5 Product description

The Invacare® Matrix® Vi and Kidabra-Vi Cushion is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The three part foam construction targets skin protection where it is needed most, and features stability, support and comfort.

**matrix**  
SEATING SERIES

**Yes, you can.®**

## 2 Safety

### 2.1 Safety information

**WARNING!**

– Do not use this product or any available optional equipment without first completely reading and understanding the user manual supplied. Invacare® product manuals are available at [www.invacare.co.uk](http://www.invacare.co.uk) or your local dealer. If you are unable to understand the warnings, cautions or instructions please contact a health care professional, dealer or technical personnel before attempting to use this equipment – otherwise, injury or damage may occur.

**WARNING!**

Invacare® products are specifically designed and manufactured for use in conjunction with Invacare® accessories. Accessories designed by other manufacturers have not been tested by Invacare® and are not recommended for use with Invacare® products.

**IMPORTANT!**

The information contained in this document is subject to change without notice.

- Check all parts for shipping damage and test before using.
- In case of damage, do not use.
- Contact Invacare® / Carrier for further instructions.

### 2.2 Symbols on the product

	Refer to user manual		Do not bleach
	CE conform		Do not iron
	Tumble dry low heat		Machine wash temperature max. 60 °C
	Do not dry clean		

## 3 Setup

### 3.1 Installing the cushion

1. Position the cushion so the material identification tags (not shown) face the rear.
2. Align the hook and loop strips (not shown) on the bottom of the cushion with those on the wheelchair seat.
3. Position the cushion onto the wheelchair as far back as possible.
4. Press down to ensure the hook and loop strips are firmly engaged.

## 4 Usage

### 4.1 Safety information

**WARNING!**

**Risk of injury**

- Installing a cushion on a wheelchair may affect the center of gravity of the wheelchair and may cause the wheelchair to become unstable, potentially resulting in injury.
- After any adjustments, repair or service and before use, make sure that all attaching component parts are secure.
- Ensure the seat cushion is properly secured to the wheelchair before using. Otherwise injury or damage may occur.
- Make sure that the locking clips are installed to ensure the safety of the user and to prevent the cushion from lifting.

**WARNING!**

**Risk of injury**

- Skin condition should be checked very frequently after the installation of any new seating system.
- Your therapist and physician should be consulted if you have any questions regarding individual limitations and needs.
- Working with your therapist, physician, and equipment supplier is the best way to ensure that a seating choice matches your individual needs.
- As the needs of the individual become more complex, the seating evaluation becomes more important.
- The cushion is designed to be used only at room temperature. If it is left at lower temperatures, always allow it to warm to room temperature before using. Sitting on extremely cold or hot surfaces can cause skin damage.
- Avoid sharp objects or exposure to excessive heat or open flame.

### 4.2 Using the cushion

1. Remove all packaging before use.
2. Place the cushion on wheelchair or chair with grey non-slip base facing down.

## 5 Maintenance

### 5.1 Inspection

**WARNING!**

– Do not continue to use this product if problems are discovered. Corrective maintenance can be performed at or arranged through your equipment supplier.

1. Visually inspect parts including drop hooks, hardware, upholstery materials, foams (if accessible), and plastics, for deformation, corrosion, breakage, wear or compression and replace if necessary.

### 5.2 Cleaning and care

**IMPORTANT!**

All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.

- For further information on decontamination in health care environments, please refer to 'The National Institute for Clinical Excellence' guidelines on Infection Control [www.nice.org.uk/CG139](http://www.nice.org.uk/CG139) and your local infection control policy.

## Cleaning the outer cover (Matrix® Vi)

1. Wipe clean the outer cover with a damp cloth.  
or
2. Carefully remove the cover by unzipping the zipper found on the rear of the cushion.
3. Close the zipper prior to laundering.
4. Machine wash at 60° C using the delicate cycle (Instructions on label).

**!** **IMPORTANT!**  
Washing at higher temperatures will cause shrinkage.

## Drying the outer cover

1. Hang cover from a line or bar and drip dry in a clean indoor environment.  
or  
Tumble dry on a low heat setting.

**!** **IMPORTANT!**  
– Tumble dry setting must not exceed 40 °C.  
– Do not tumble dry for longer than 10 minutes.  
– Dry thoroughly before re-fitting to the foam.

## Cleaning the inner cover

**!** **IMPORTANT!**  
**Risk of damage to the product**  
– Do not remove the inner cover from the foam core.

1. Wipe clean the inner cover with a damp cloth.

## Replacing covers

1. Carefully remove the cover by unzipping the zipper found on the rear of the cushion.
2. Place new cover onto the foam core.

**!** **IMPORTANT!**  
– Ensure that the corners of the foam core are positioned correctly into the corners of the cover.

## 6 After Use

### 6.1 Storage

**!** **IMPORTANT!**  
– Store cushions in a dry environment.  
– Store cushions within a protective cover.  
– Store items on clean, dry, off-flooring free from sharp edges to avoid any possible damage.  
– Never store other items on top of a cushion.  
– Do not store cushions next to radiators or other heating devices.  
– Protect cushions from direct sunlight.

1. Refer to the storage and shipping conditions in section 7.3 Environmental Parameters, page 3.

### 6.2 Re-use

The product is suitable for repeated use. The number of times it can be used depends on how often and in which way the product is used.

1. Before reuse, clean the product thoroughly, see section 5.2 Cleaning and care, page 2.

### 6.3 Disposal

The disposal and recycling of used devices and packaging must comply with the applicable local legal regulations.

1. Ensure that the cushion is cleaned prior to disposal to avoid any risk of contamination.

## 7 Technical data

### 7.1 General data

Product	Kidabra-Vi cushion	Matrix® Vi cushion	Matrix® Vi heavy duty cushion
Total height [mm]	38	57	57
Total width [mm]	300 – 360	360 – 510	460 – 760
Total depth [mm]	300, 360	380 – 510	460 – 660
Maximum user weight [kg]	68	136	272
Weight complete cushion [kg] <sup>1)</sup>	0.5 – 0.7	0.9 – 2.0	1.6 – 4.0
Weight outer cover [kg] <sup>1)</sup>	0.19	0.25 – 0.44	0.4 – 1.0

<sup>1)</sup> Weights can vary depending on size of cushion ordered.

Size charts and local requirements / certification available upon request from Invacare®.

### 7.2 Materials

Foam top layer	High resilient Polyurethane foam
Foam bottom layer	High resilient Polyurethane foam
Cover	Knitted polyamide, coated with an anti-microbial polyurethane coating. Substrate treated with a fluorocarbon to give water/stain resistance.
All cushion components are Latex-free.	

### 7.3 Environmental Parameters

Operating conditions	
Ambient temperature	10 - 35 °C
Relative humidity	30% - 75% Non-Condensing
Atmospheric pressure	70 - 106 kPa
Storage and shipping conditions	
Ambient temperature	-40 - 70 °C
Relative humidity	10% - 100%
Atmospheric pressure	50 - 106 kPa

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Invacare UK Operations Limited