





SEAT HEIGHTS ACTION³ JUNIOR



Seat Angle 5°	Rear Support		Headset Support	
	20"	22"	Position 1 (high)	Position 2 (low)
41	4	5	X	3
43.5	3	4	2	X
46	2	3	1	2



Seat Angle 1.5°	Rear Support		Headset Support	
	20"	22"	Position 1 (high)	Position 2 (low)
41	3	4	X	3
43.5	2	3	2	X
46	1	2	1	2



Seat Angle 8.5°	Rear Support		Headset Support	
	20"	22"	Position 1 (high)	Position 2 (low)
41	5	X	X	3
43.5	4	5	2	X
46	3	4	1	2

Delivered configuration
 Not technically possible

