



Product Service

# CERTIFICATE

No. Z1 07 10 12941 105

**Holder of Certificate:** Invacare Deutschland GmbHKleiststrasse 49  
32457 Porta Westfalica  
GERMANY**Certification Mark:****Product:**

Electrical wheelchairs class B

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. The certification mark must not be altered in any way. See also notes overleaf.

**Test report no.:**

71326670

**Date,** 2007-10-25

Page 1 of 2



Product Service

**CERTIFICATE**

No. Z1 07 10 12941 105

**Model(s):** **Electrically driven wheelchair**  
**Type: Typhoon II**

**Parameters:**

model:	Typhoon II, middle-wheel drive
class of use:	B
total height:	1070 mm – 1370 mm
total length:	1095 mm
total width:	650 mm - 695 mm
stroke:	300 mm
obstacle height:	80 mm
max. speed:	6 km/h / 10 km/h 12 km/h
max. load:	150 kg
max. safe slope:	11,5° (20%, in the standard seat-position, as specified of the manufacturer)
driving motor:	Invacare Deutschland GmbH, construction: permanently excited multipolar three-phase current motor, article-no. 1433561 & 1433562
drive of the adjustment of the seat inclination:	REAC AB, type: RE23
power electronics:	Dynamic Controls Ltd., type: DX-GB-ASG
batteries:	MK Battery, 2x(12 V, 60 Ah (20hr)), type: M34 SLD G
battery charger:	CTE Corporation, model: BAT-GCO812, 24 V DC, 8 A

**Tested according to:** DIN EN 12184:1999

**Production Facility(ies):** 12941