



Product Service

CERTIFICATE

No. Z1 05 07 12941 099

Holder of Certificate: Invacare Deutschland GmbH

Kleiststrasse 49
32457 Porta Westfalica
GERMANY

Certification Mark:



Product:

**Wheelchair for handicapped
Electrically stand-up wheelchair**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. See also notes overleaf.

Test report no.:

70102920



186734

Date, 2005-07-29

Page 1 of 2



Product Service

CERTIFICATE

No. Z1 05 07 12941 099

Model(s): **Electrically stand-up wheelchair**
Type: Dragon Vertic

Parameters:

class of use:	class B, rear-wheel drive	
total height:	1155 mm - 1680 mm	
total length:	1115 mm	
total width:	625 mm	
max. load:	100 kg	
electro-motor:	Invacare, type: SSG0220	
electronics:	Dynamic Controls Ltd., type: Shark (remote: DK-REMD21, power-module: DK-PMB21)	
batteries:	MK Battery, 2x (12 V, 60 Ah (20hr)), type: M34SLDG	
battery charger:	CTE Corporation, model: BAT-GCO812, 24 V DC, 8 A	
	sitting position	stand-up position
max. safe slope:	10°	3°
max. obstacle height:	40 mm	0 mm
max. speed:	6 km/h	1,2 km/h

Tested according to: EN 12184:1999
 ISO/CD 7176-20:2003

Production Facility(ies): 12941