

**Intertek****ETL SEMKO****EMC Pre-Compliance Summary Sheet**

<b>Project No.:</b> 2003-07-579	<b>Date:</b> 23 July 2003
<b>Company:</b> Invacare UK Limited	
<b>Product:</b> Mirage (rabbit) Power Chair, Shark Controller, 3.6 Amp Huafeng Motors	
<b>Standard:</b> EN 12184 : 1999	
<b>Test</b>	<b>Comments</b>
Radiated Emissions	Refer to plots.
Radiated Immunity 10V/m, 80% AM, 80-1000 MHz	Vertical Polarisation - No malfunctions observed. Horizontal Polarisation - No malfunctions observed.
ESD ± 6 kV Contact ± 8 kV Air	No malfunctions observed

Test Engineer: Ginnie Crute



**Note:** The information on this sheet in no way constitutes a report and no details of the equipment under test or test results are kept in Interteks records.

ITS Testing &amp; Certification Ltd

Blackwood EMC Facility, Unit 8 Woodfieldside Business Park,  
Pontllanfraith, Blackwood, Gwent, NP12 2DGTel: + 44 (0)1495 229219 Fax: +44 (0)1495 228331 Web: [www.uk.intertek-etlsemko.com](http://www.uk.intertek-etlsemko.com)

Registered No 3272281 Registered Office: 25 Savile Row London W1S 2ES

