



Product Service

# CERTIFICATE

No. Z1 09 05 11048 011

**Holder of Certificate:** Ulrich Alber GmbHVor dem Weißen Stein 21  
72461 Albstadt  
GERMANY**Production Facility(ies):** 11048**Certification Mark:****Product:** Mobilisation aids**Model(s):** Mobile attendant-operated stair-climbing device type: scalacombi® eco S34

<b>Parameters:</b>	max. user weight:	120 kg
	total width:	520 mm
	total height:	1360 mm
	rated voltage:	24 V DC
	drive power:	200 W
	degree of protection:	IPX3
	speed:	12 steps per minute
	theoretical number of steps:	200
	batteries:	2 x 12VDC, 3.4 Ah
	battery charger:	BC 24005LA, 90-265 VAC, 24 VDC, 0.5A

**Tested according to:** ISO 7176-23:2002  
ISO 7176-14:2008  
DIN 12182:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

**Test report no.:** 71348307**Date,** 2009-05-12

Page 1 of 1

