



Trusted Solutions,
Passionate People

Power Wheelchair
Ireland Prescription Form
PRONTOM41CAB010117

Tel: (00353) 1 8107084 Fax: (00353) 1 8107085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JANUARY 2017

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



Pronto M41 Ireland

Cab Seating

Max user weight 136kg (21.4 stone)

Standard Features



- Black Cab Seat
- Seat width 406 mm (16") or 457 mm (18")
- Backrest height 480 mm (19")
- Flip up armrests
- Padded height adjustable headrest
- Manual setttable reclining backrest recliner by lever
- Adjustable flip-up centre footboard
- Standard Remote MK5 SPJ joystick with standard bracket
- Candy Red shroud
- Solid Insert/Puncture Proof Grey Mid Wheel (255mm x 80 mm, 10" x 3")
- Disengage high torque motors
- Kerb Climbing range 0-50 mm
- Turning circle 1160 mm
- Battery 32Ah AGM pair
- Metal Buckle Lap Belt (Adjustable from Both Sides)



* Configurations may vary from the image shown

Standard Base Price €2,922

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code		Description	IRELAND	
AQB0000		Pronto M41 Basic Model	◆	
AQZ0000		Pronto M41 Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM	POA	
SEATING TYPE				
AQB1172	X	Cab Seat (black)		€ 2,922
SEAT SIZES				
AQB0341		Seat width 406 mm (16"), with Seat Depth 406 mm (16")		○
AQB0346		Seat width 457 mm (18"), with Seat Depth 457 mm (18")		○
BACKREST HEIGHT				
AQB0448	X	480mm (19")		◆
ARMREST				
AQB0609	X	Flip-up armrests, height and width adjustable		◆
SEAT ACCESSORIES				
HEADREST				
AQB1491	X	Headrest		◆

SEAT ADJUSTMENTS			
BACKREST RECLINER			
AQB0959	<input checked="" type="checkbox"/>	Man. settable reclining angle by lever (95°-115°)	◆
SEAT TO FLOOR HEIGHT			
AQB0364		Seat height: 470 mm (18.5")	○
AQB0322		Seat height: 495 mm (19.5")	○
AQB0321		Seat height: 520 mm (20.5")	○
AQB0323		Seat height: 545 mm (21.5")	○
AQB0359		Seat height: 570 mm (22.5")	○
CENTRE MOUNTED LEGREST			
AQB0595	<input checked="" type="checkbox"/>	Adjustable flip-up footboard (no hangers)	◆
REMOTES			
AQB1840	<input checked="" type="checkbox"/>	Standard Remote MK5 SPJ joystick	◆
AQB1207		Standard Bracket RH	○
AQB1208		Standard Bracket LH	○
DRIVE UNIT OPTIONS			
SHROUD COVER COLOUR			
AQB0227	<input checked="" type="checkbox"/>	Candy Red	◆
WHEELS			
AQB0781	<input checked="" type="checkbox"/>	Centre Wheel 255mm x 80mm (10"-3") Solid Insert/Puncture Proof Grey	◆
PERFORMANCE OPTIONS			
AQB1282		Motor type 6km/h (4mph)	◆
AQB1294		Motor type 3km/h (2mph) Indoor Programming	○
AQB1261	<input checked="" type="checkbox"/>	Battery 32Ah AGM pair	◆
AQB1343	<input checked="" type="checkbox"/>	Charger 8A (200-250V)	◆
SAFETY OPTIONS			
AQB1329	<input checked="" type="checkbox"/>	Metal Buckle Lap Belt (Adjustable from Both Sides)	◆
TOOLKIT			
AQB1333	<input checked="" type="checkbox"/>	Toolkit (incl bit holder, bits and fork spanners)	◆
			Total Price



Powerchair Alternative Options

Need a different option not on the form?

Take a look at our Powerchair Alternative Options form, available as standard.

Powerchair Alternative Option required:



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):