

# Technical Information

## - Pressure Area Care Products

Product	Warranty	Fire Testing	Grade Ref & Colour
<i>Invacare</i> <b>Softform®</b>			
<b>Premier Active</b>	4 Years cover 8 Years foam 2 Years insert 2 Years pump	BS 5852: Crib 5	RX39/120 Pink RX39/200 Blue
<b>Premier MaxiGlide</b>	4 years cover 8 years foam	BS 5852: Crib 5 BS 7175: Crib 5	GR. RX 39/120 Pink GR. RX 39-200 Blue
<b>Premier</b>	10 years foam 5 years cover	BS 7177: Crib 7	Gr. 39/120 Pink Gr. RX 39/200 Blue
<b>Premier Visco</b>	4 Years *	BS 5852: Crib 5	Gr. RX 39/120 Pink Gr. RX 39/200 Blue Gr. VC5580 White
<b>Excel</b>	2 Years	BS 7177: Crib 5	Gr. RX 39/120 Pink Gr. RX 39/200 Blue
<b>Heel Pad</b>	1 Year	BS 5852: Crib 5	Gr. RX 33/125 White
<b>Odstock Wedge</b>	1 Year	BS 5852: Crib 5	Gr. RX 33/125 White
<i>Invacare</i> <b>Propad®</b>			
<b>Premier Overlay</b>	3 Years**	BS 7177: Crib 5	Gr. RX36/125 Beige
<b>Original Overlay</b>	3 Years**	BS 7177: Crib 5	Gr. RX36/125 Beige
<b>Visco Overlay</b>	1 Year	BS 5852: Crib 5	Gr. VC5580 White Gr. RX36/125 Beige
<b>Premier Cushion</b>	1 Year	BS 5852: Crib 5 BS 7175: Crib 5	Gr. RX36/125 Beige Gr. RX39/200 Blue
<b>Original Cushion</b>	1 Year	BS 5852: Crib 5 BS 7175: Crib 5	Gr. RX36/125 Beige
<b>Low Profile Cushion</b>	1 Year	BS 5852: Crib 5 BS 7175: Crib 5	Gr. RX36/125 Beige
<b>Visco Cushion</b>	1 Year	BS 5852: Crib 5 BS 7175: Crib 5	VC5580 RX39/200 Blue
<b>Leg Trough</b>	1 Year	BS 5852: Crib 5	RX39/200 Blue
<i>Invacare</i> <b>Flo-tech®</b>			
<b>Flo-tech® Cushions</b>	3 Years**	BS 7176: 1995 Crib 5	Gr. FG 70/180 Grey
<b>Flo-tech® Visco Cushion: Soft</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. LM 85/115 Blue
<b>Flo-tech® Visco Cushion: Medium</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. LM 85/150 Blue
<b>Flo-tech® Visco Cushion: Firm</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. LM 85/190 Blue
<b>Flo-tech® Visco Cushion: Very Firm</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. LM 85/250 Blue
<b>Flo-tech® Back Support</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. FG 70/180 Grey
<b>Flo-tech® Back Support Lumbar</b>	3 Years**	BS 5852: 1990 Crib 5	Gr. FG 70/180 Grey VC55085 White
<b>Lumbar Pad</b>	1 Year	BS 5852: 1990 Crib 5	Gr. RX 39/200 Blue VC55085 White



**Yes, you can.®**