

Invacare® | **XLT Max**




### Stability, Rigidity, Security: Be Active!

The new *Invacare XLT Max* is a bariatric, rigid active wheelchair ideal for day-to-day use. Built with a titanium frame, it is designed to support a user weight of max. 180 kg with seat width possibilities of up to 600 mm. The **XLT Max** offers individuals various levels of support enabling greater mobility and a more active lifestyle.



**Yes, you can.®**



## The Perfect Fit

Measured and constructed to fit the requirements of each person, the **Invacare XLT Max** is an ideal chair for those previously unable to take advantage of an active wheelchair. Built for maximum security and support in combination with good driving performance and manoeuvrability, the **XLT Max** is the only active rigid bariatric wheelchair on the market today.

## Simply Compatible

The **XLT Max** is a rigid titanium wheelchair offering stability and the utmost rigidity. The well-designed fit of the wheelchair helps people with a weight of up to 180 kg retain and improve mobility and independence in their daily lives. Individuals will benefit from the extensive and readily available array of Invacare options which can easily be fitted to the **XLT Max**.

## Features and Options



### Invacare Flo-tech Contour

Ideal bariatric cushion solution. Improved posture and increased support surface area significantly contributes to preventing pressure ulcers.



### Footplate converter

Divided into three parts, the footplate comes with a removable mid section, to facilitate transfers.



### Anti-tipper device

The anti-tipper device is height adjustable.

## Technical data



**Invacare XLT Max** 505 / 555 /  
605 mm



500 / 550 /  
600 mm



440 - 520 mm



400 - 490 mm



380 - 480 mm  
270 - 370 mm high mounted



0 - 14°



**Invacare XLT Max** ± 10°



SW + 210 mm



790 - 1100 mm



930 - 1080 mm



15.5 kg



180 kg



**Invacare XLT Max** 11.1 kg



Crash-tested  
ISO 7176-19



The **Invacare XLT Max** complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: [www.invacare-re.com](http://www.invacare-re.com) under Rea wheelchairs/Technical downloads.

INVACARE Limited  
Pencoed Technology Park - Pencoed -  
Bridgend CF35 5AQ  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)

INVACARE Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2011 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare reserves the right to alter product specifications without prior consultation.

XLT Max - UK - 07/2011 - 1538967



**Yes, you can.®**



Blue



Red



Black



Grey



Blasted Titanium



Orange