

Invacare® | **Commode chairs**



**160kg
maximum user
weight**

Supplied with leg rests as standard



2 Year Warranty

Invacare **Glideabout commode**

- Stable, robust frame made from powder-coated steel
- High corrosion protection
- Space-saving push bar for optimal positioning over the toilet
- Swing-away armrests and swing-away, detachable footrests
- Stable commode seat with integrated toilet pan holder
- Supplied with 4 braked castors as standard



Yes, you can.®



Invacare Omega height adjustable

- Height adjustable commode chair
- Space-saving: can be stacked
- Seat height adjustable in 7 steps from 510 - 660 mm
- Attractive design: fits well into a domestic setting
- Easy to clean






Invacare Omega Piante

- Foldable, Space saving commode
- Ideal for infrequent or travel use
- Contoured backrest
- Complete with sanitary pan

Technical data

Order numbers	Invacare <i>Glideabout commode</i>	Invacare <i>H450LA height adjustable</i>	Invacare <i>Omega H407</i>
Commode chair, mobile	1525771		
Omega mobile commode chair with castors			
Static commode		1471166	
Omega height adjustable commode chair			1471010 (2 pack)
Dimensions			
Seat width	450 mm (A)	495 mm (A)	440 mm (A)
Seat depth	450 mm (B)	455 mm (B)	420 mm (B)
Seat height	480 mm (C)	510 - 660 mm (C)	450 mm (C)
Total width	575 mm (D)	670 mm (D)	580 mm (D)
Total depth	810 mm (E)	630 mm (E)	560 mm (E)
Total height	945 mm (F)	900 - 1050 mm (F)	910 mm (F)
Depth without footrests	545 mm (G)	–	–
Footplate depth	155 mm (H)	–	–
Height to armrests	755 mm (I)	700 - 850 mm (H)	660 mm (H)
Height, footrests to seat	345 - 460 mm (J)	–	–
Height, ground to lower seat edge: with holder for toilet pan	420 mm (K)	380 - 530 mm (I)	–
without holder for toilet pan	–	440 - 590 mm (J)	–
Height, floor to cover plate	525 mm (L)	535 - 685 mm (K)	475 mm (K)
Width between the armrests	470 mm (M)	510 mm (G)	460 mm (G)
Width, cut out in seat plate	220 mm	215 mm	215 mm
Depth, cut out in seat plate	260 mm	250 mm	255 mm
Ø Castors	5" (127 mm)	–	–
More info			
Weight	14,0 kg	8,8 kg	8,5 kg
Load capacity	160kg	130 kg (21 stone)	80 kg (13.5 stone)
Colour	Black anthracite	Lavender blue / Blue	Lavender blue
Warranty	2 years	2 years	2 years

 Invacare *Glideabout commode*

 Invacare *Omega Height adjustable*

 Invacare *Omega Pilante*




Invacare Styxo 9630E

The *Invacare 9630 Styxo* combines two functions in one product. It can either be used as toilet frame over the toilet or as a toilet commode with usage of the toilet pan. The toilet pan with lid and the splash guard are included in the delivery.

Product Description:

- Easy to handle seating height adjustment in 7 steps from 410mm to 560mm
- Complete with elongated toilet seat and lid, pan with lid and splash guard
- Tool less assembly and easy to clean
- High load capacity: 160kg (25 Stone)




Invacare Adeo C407A

- Height adjustable 640 - 740 mm
- Fits most standard toilets
- Complete with armrests

Technical data

Order numbers	<i>Invacare Styxo</i>	<i>Invacare Adeo</i>
Toilet frame	9630E	C407A
Dimensions		
Seat width	360 mm	590 mm
Seat depth	420 mm	440 mm
Seat height	410-560 mm	640 - 740 mm
Total width	620 mm	590 mm
Total depth	490 mm	440 mm
Total height	600 - 750mm	740 mm
More info		
Weight	5,5 kg	2kg
Load capacity	160kg	100kg (16 stone)
Colour	grey	Chrome
Warranty	2 years	2 years

Invacare Styxo 9630E



Invacare Adeo C407A

