



küschall **R33**

Active wheelchair

New!

Precision is essential for *küschall*, and pushes us to bring more and more information to you. This is the reason why this new portfolio edition includes the weight of each part.

Thanks to these details, you will be able to compare the weight of all options and make a safe choice; you will then easily estimate the whole weight of the *R33* wheelchair you have set.

Except any other indication, all weights indicated are based on seat width 40cm and seat depth 40cm. Options weights are calculated by pair when requested (e.g. for the wheels, push handles, brakes...) and with their related fixations.

R33

Contents

Frame.....	4
Backrest.....	6
Seating	8
Footrest.....	9
Side guards.....	10
Cushions	11
Colours	12
Front wheels	13
Rear wheels	14
Brakes	17
Safety options.....	18
Accessories	19
<i>Küschall</i> lifestyle collection.....	20

Frame

Frame angle

DDR0237

Frame angle 75° to the ground



DDR0238

Frame angle 90° to the ground



Seating possibility

DDR0626

Low Impact System



DDR0627

Fixed seating



Rear seat-to-floor height (SHh)

With Low Impact System



Spring, user weight <70kg			Spring, user weight 70 - 85kg			Spring, user weight 85 - 100kg		
SHh43	SHh45	SHh47	SHh43	SHh45	SHh47	SHh43	SHh45	SHh47

With Fixed Seating

SHh40	SHh41	SHh42	SHh42	SHh43	SHh44	SHh44	SHh45	SHh46	SHh46	SHh47	SHh48

Backrest

Backrest material

DDR0730	Aluminium	
	<i>Short: 320 grams</i> <i>Large: 380 grams</i>	
DDR0731	Titanium	
	<i>Short: 230 grams</i> <i>Large: 240 grams</i>	

Push handles

DDR0750	Without push handle	
DDR0702	Push handles fixed, long	
	<i>240 grams the pair (without back tube)</i>	
DDR0751	Push handles fixed, short	
	<i>160 grams the pair (without back tube)</i>	
DDR0704	Push handles, foldable	
	<i>1'120 grams the pair (incl. tubes and fixations)</i>	

Backrest upholstery

DDR0701	Black backrest, tension adjustable with <i>Velcro</i> [®]	
	<i>560 grams (including Velcro[®] straps)</i>	
DDR0735	Aero backrest, tension adjustable with <i>Velcro</i> [®]	
	<i>600 grams (including Velcro[®] straps)</i>	
DDR0748	Alcantara black backrest, tension adjustable with <i>Velcro</i> [®]	
	<i>681 grams (including Velcro[®] straps)</i>	
DDR0746	Alcantara backrest, tension adjustable with <i>Velcro</i> [®] special edition red & dark grey	
	<i>681 grams (including Velcro[®] straps)</i>	
DDR0747	Alcantara backrest, tension adjustable with <i>Velcro</i> [®] special edition dark grey & light grey	
	<i>681 grams (including Velcro[®] straps)</i>	
DDR0749	Alcantara backrest, tension adjustable with <i>Velcro</i> [®] special edition red & black	
	<i>681 grams (including Velcro[®] straps)</i>	
DDR0752	Alcantara backrest, tension adjustable with <i>Velcro</i> [®] special edition black & dark grey	
	<i>681 grams (including Velcro[®] straps)</i>	


Seating

Seat upholstery

DDR0610

Standard seat upholstery


220 grams (including fixations)



DDR0611



Tension adjustable

560 grams (including Velcro® straps and fixations)





Footrest

Footrest type

DDR1510	Standard footrest tube, rigid	
	<i>370 grams</i>	
DDR1570	Carbon footplate, angle adjustable	
	<i>790 grams (including fixations)</i>	

Footrest tube material

DDR1508	Aluminium	
	<i>360 grams</i>	
DDR1509	Titanium	
	<i>210 grams</i>	

Side guards

Clothes guards / Mud guards

DDR1401 Carbon clothes-guards, flip back

*Short: 420 grams the pair (incl. fixations)
Large: 540 grams the pair (incl. fixations)*



DDR1408 Carbon mud guards, removable
(This option adds 2 cm to the total width of the wheelchair)

*Short: 660 grams the pair (incl. fixations)
Large: 824 grams the pair (incl. fixations)*









Siderests

DDR1416 Siderests height adjustable (with tools),
DDR1421 removable, short or long

*Short: 1'580 grams the pair
Long: 1'700 grams the pair*



Cushions






DDR600	Seat cushion	NO PICTURE AVAILABLE
	<i>510 grams (SB40 cm*ST40 cm)</i>	
DDR0660	Flo-tech Lite cushion	
	<i>1'020 grams (SB40 cm*ST42 cm)</i>	
DDR0603	Alcantara cushion, special edition red & dark grey	
	<i>460 grams (SB40 cm*ST36 cm)</i>	
DDR0604	Alcantara cushion, special edition dark grey & light grey	
DDR0607	Alcantara cushion, special edition black & red	
DDR0608	Alcantara cushion, special edition black & dark grey	
DDA0605	Alcantara cushion, black	
	<i>460 grams (SB40 cm*ST36 cm)</i>	

Colours

Frame colours




Glossy black		Black matt	
Royal blue		Royal red	
Purple fire		Silver grey	
Emerald green		Midnight blue	
Umber glaze		Ivory white	
Smokey grey		Night shimmer	
Mirror polished			

CPS-colours

Techno-anthracite		Lemon	
B-blue		Comet red	
Artic blue			

The Colour Power Set (CPS) applies on the following parts and options: Starec & High Performance rear wheel hubs with quick release protection, 24" Starec rear wheel rims, Starec front wheels, forks, castor housing caps, front seat module caps, backrest adjustment plates, L.I.S.[®] centre part.

Front wheels

DDR0800	Supersport	
	<i>3": 212 grams the pair</i>	
	Starec wheel, in various CPS colours	
	<i>4": 364 grams the pair</i> <i>5": 536 grams the pair</i>	
DDR0863	Comfort wheel with black rubber	
	<i>5": 780 grams the pair</i>	

Rear wheels

Axle type

DDR0278

Carbon Axle

230 grams



Rear wheel type

DDR0915

kuschall Starec wheel

24": 1'592 grams the pair



DDR0908
DDR0918

High Performance wheel
(standard hub is black – if no CPS choice)

24": 1'356 grams the pair
25": 1'372 grams the pair



DDR0913

Mountain Bike wheel

24": 1'712 grams the pair



DDR0917

Spinergy wheel

24": 1'476 grams the pair
25": 1'508 grams the pair





Rear wheels tyres

<p>DDR1001 Light wheel 7 bar</p>  <p>24": 816 grams the pair</p>	<p>DDR1002 Profile tyre 7 bar</p>  <p>24": 624 grams the pair</p>
<p>DDR1007 Solid, light wheel</p> <p style="color: red; text-align: center;">NO PICTURE AVAILABLE</p> <p>24": 1'032 grams the pair</p>	<p>DDR1004 Solid tyre, profile</p>  <p>24": 1'140 grams the pair</p>
<p>DDR1006 Solid tyre (KIK Type), black</p>  <p>24": 904 grams the pair</p>	<p>DDR1010 Marathon Evolution Plus</p>  <p>24": 1'024 grams the pair 25": 1'076 grams the pair</p>
<p>DDR1011 Schwalbe speedrun</p>  <p>24": 356 grams the pair 25": 996 grams the pair</p>	

Handrims

<p>DDA1102 Supergripp</p>  <p>24": 1'336 grams the pair</p>	<p>DDA1104 Steel</p>  <p>24": 1'276 grams the pair</p>
<p>DDA1106 Titanium</p>  <p>24": 808 grams the pair</p>	<p>DDA1113, DDA1114 Light Alu</p>  <p>24": 716 grams the pair 25": 744 grams the pair</p>


Spoke guards

DDR1200	Transparent	
	<i>24": 660 grams the pair</i>	
DDR1207	«Original» Carbon	
	<i>24": 720 grams the pair</i>	
DDR1202	küsch	
	<i>24": 680 grams the pair</i>	
DDR1209	Special R33	
	<i>24": 720 grams the pair</i>	

Brakes

DDR1300 Push brakes, levers bent

380 grams the pair



DDR1302 Performance brakes

420 grams the pair





DDR1303 Active brakes

300 grams the pair











Safety options

DDR1702 DDR1703 DDR1704	Antitipper, height adjustable, active and inactive position, removable, right, left or both sides <i>Short: 670 grams (including fixations)</i> <i>Long: 700 grams (including fixations)</i>	
DDR1730	Seat belt <i>230 grams</i>	
DDR1710 DDR1711	Passive lights, assembled or enclosed <i>60 grams the pair</i>	NO PICTURE AVAILABLE

Accessories

DDR1728	Wheel bag for 24", 25" and 26" rear wheels	
	<i>750 grams</i>	
DDR1715	Air pump	
	<i>100 grams</i>	
DDR1716	AirMan	
DDR1717	Tool kit	
	<i>60 grams</i>	
DDR1750	Axle protection	NO PICTURE AVAILABLE
DDR1739	Front frame protectors (2 pieces)	
	<i>80 grams the pair</i>	

Kuschall lifestyle collection

1498701 1498702 1498734 1498735 1498776	T-Shirt R33 men or women (available in S, M, L, XL, XXL)	
1498736/1498778 1498737/1498779 1498738/1498780 1498739/1498780 1498777	Polo R33 men (available in S, M, L, XL, XXL) or women (available in S, M, L, XL)	
1531382 / 1531383 1531463 / 1531464 1531189 / 1531379 1531380 / 1531381	<i>Kuschall</i> Black T-Shirt men (available in size M to XXL) or women (available in size S to XL)	
1513357 1513358 1513359 1513360 1513361	<i>Kuschall</i> Jacket (available in S, M, L, XL, XXL)	
1518297 1518298 1518299 1518300 1518301	<i>Kuschall</i> Sweat-shirt (available in S, M, L, XL, XXL)	
1531187	<i>Kuschall</i> Baseball Cap Black & Silver	
1486578	<i>Kuschall</i> Pen	
1486577	<i>Kuschall</i> Carbon Pen Set	

1486614 *Kuschall* Umbrella



1532047 *Kuschall* Bag-pack



1486620 *Kuschall* measuring tape blue



