Invacare True Track Plus





5 year warranty



Go further and faster in total control

The new generation *True Track Plus* motor gives more driving range and power for clients who want the best performance in speed and control. *True Track Plus* immediately responds to the client's commands and provides more freedom in driving by easily and smoothly overcoming obstacles. Clients will appreciate the speed and discretion that comes with this silent motor.





Unrivalled efficiency, power and torque

Improved efficiency gives more range and power, with easier slope climbing at slower speeds with up to 40% more torque!

Power and immediate response are maximised for clients of all levels of experience, with many options of driver controls so more experienced clients can enjoy a faster and more dynamic drive.

Absolutely silent and maintenance free

No friction, no wearing parts, no carbon brushes – no maintenance.

This option completely removes the noise and vibrations that occur on conventional geared motors. The gearless, brushless motor provides maximum discretion for the customer.



Speed choice

Adapt to your client's needs and easily programme all speed choices without loss of torque.



Smarter design

The latest in engineering means the **True Track Plus** motors are high performance and maintenance free.



ACS2 controls

ACS2 controls offer a number of input options.

True Track Plus is available for Invacare Storm®4 X-plore.

.....

5 year warranty

Technical data

True Track

Plus



600 W



6 / 10 / 13 km/h



59 km⁻¹

1/ according to ISO 7176-4 on Storm4 X-plore

Invacare Limited
Pencoed Technology Park - Pencoed Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.com

Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS - County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com

Email: Ireland@invacare.com www.invacare.ie

© 2014 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

True Track Plus - UK - 09/2014 - 1579018

