

# Invacare® | **Zimma™ Frames**



## **Trustworthy, stable, and reliable**

The NEW *Invacare Zimma Frame* is a robust fixed walking frame designed for use in both domestic and hospital environments.

Both types are available in a range of sizes to meet the needs of different individuals. The frames are also fully height adjustable making them adaptable to different situations.

Available with or without wheeled extension legs, the *Invacare Zimma Frame* is trustworthy, stable, and reliable.



**Yes, you can.®**



**Features and Options**

**Range of Types and Sizes**

The *Invacare Zimma Frame* is available in both Domestic and Hospital frame types. Both types are available in small, medium and large sizes to suit a range of different patient requirements.

All *Invacare Zimma Frames* are height adjustable by way of a simple push button giving a series of fixed settings.

**Wheeled Extension Legs**

The wheeled extension legs are available for all *Invacare Zimma Frames*.

They are simple to fit, giving a quick and easy way of converting the frame to a wheeled version.



Range of Types and Sizes



Wheeled Extension Legs

**Technical data**

<i>Invacare Zimma Frames</i>	Height (mm)	Front Width (mm)	Back Width (mm)	Depth (mm)	Product Weight	Capacity
Hospital Small	665-741	480-500	650-665	530-545	1.6 kg	160 kg
Hospital Medium	765-841	480-500	660-675	540-555	1.8 kg	160 kg
Hospital Large	865-941	480-500	670-685	550-565	2 kg	160 kg
Domestic Small	700-776	440-460	630-655	470-495	1.8 kg	160 kg
Domestic Medium	800-876	435-450	625-650	460-485	1.9 kg	160 kg
Domestic Large	900 - 976	430-445	620-645	450-475	2 kg	160 kg

Invacare Limited  
 Pencoed Technology Park - Pencoed -  
 Bridgend CF35 5AQ - United Kingdom  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County  
 Dublin - Ireland  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2011 Invacare International Sàrl  
 All rights reserved.

All information quoted is believed to be correct at time of print.  
 Invacare reserves the right to alter product specifications without prior consultation.

Invacare Zimma Frames - UK - 11/2011  
 1552285



**Yes, you can.®**

**Product Code**

**1536643 Hospital Small**  
**1536695 Hospital Medium**  
**1536696 Hospital Large**

**1536697 Domestic Small**  
**1536698 Domestic Medium**  
**1536699 Domestic Large**

**1536700 Pair of Wheeled Extension Legs**

**1552962 Hospital Small with wheels\***  
**1552963 Hospital Medium with wheels\***  
**1552964 Hospital Large with wheels\***

**1552965 Domestic Small with wheels\***  
**1552966 Domestic Medium with wheels\***  
**1552967 Domestic Large with wheels\***

\*Wheeled extension legs not fitted

