

DOLOMITE | JAZZ²

NEW



JAZZ

think
independence
www.thinkindependence.com

An innovation in independent living

The **DOLOMITE JAZZ²** is an innovative lightweight rollator designed to offer excellent support and stability, both during indoor and outdoor walks. Its ergonomic design allows a natural walking position within the frame.

This new stylish rollator includes various smart features making it the ideal walking companion. The latest member of the **DOLOMITE** family combines style with exceptional functionality.



Innovative folding frame that stands alone when folded


DOLOMITE


Yes, you can.®



Features

Smart rollator

DOLOMITE JAZZ² introduces an innovative side folding frame. The easy to operate design provides extra walking space and can be stored upright once folded. Its push handles are easily height adjustable with clear positioning marks. The memory function allows them to be lowered during transport but put back into the preset position without any adjustment. Push handles also feature an anatomic design providing pressure reduction.

Stylish

Modern, anthracite frame design combined with light grey accents, the DOLOMITE JAZZ² is a fresh and modern looking rollator which integrates easily into day to day activities. The braking cables are fully integrated into the frame providing a clean look and also prevent them from catching on objects whilst in use.



Basket and padded seat

Easily removable basket which can be used as a shopping bag. Rigid soft padded seat for great comfort and safety. (Standard feature)



Kerb climber

Integrated aid to overcome kerbs. Can be positioned left or right. (Standard feature)



Back strap

Comfortable and safe support to lean against while sitting. Available with (top picture) or without extra padding (lower picture). (unpadded backstrap supplied as standard).

Padded Back strap: **SP1551654**

Options



Large basket

Provides extra storing space with lid. Available in grey or red. (Optional)

Basket Grey: **SP1560336**

Basket Red: **SP1564104**



Cane holder

Allows safe transport of a cane or a crutch and makes it easily accessible. Easy to fit. (Optional)

SP1512912



Tray

Removable tray to safely transport items. (Optional)

SP1546264

Technical data

JAZZ ² 610	598 mm	689 mm	510 mm	750 – 1000 mm	610 mm	7,8 kg
JAZZ ² 510	598 mm	689 mm	510 mm	680 – 800 mm	510 mm	7,5 kg
JAZZ ² 610 / 510	454 mm	200 mm	150 kg (24 stone)	215 mm		

Frame colour



Anthracite grey

Product Codes: (with unpadded backstrap included)

JAZZ² 610 - 1559202-46-746

JAZZ² 510 - 1559201-46-746

Including basket and unpadded backstrap

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

Jazz² - UK - 01/2013 - 1558677



Yes, you can.®

