



INCREASED CLINICAL FUNCTIONALITY

Different people have different clinical seating requirements - to meet the individual needs we have developed the Invacare® **Flex3** system as a supportive and flexible seat and backrest system.

The **Flex3** system offers combinations of backrest types and backrest cushions from little support to highly ergonomic positioning for people that require low as well as a high level of comfort and support.

COMPATIBILITY BETWEEN PRODUCT LINES

With the same backrest attachment on all wheelchairs, the **Flex3** system offers the opportunity to combine backrest cushions and backrest types across the range of Passive wheelchairs. In addition it allows backward compatibility and enhances refurbishment.

INCREASED CLINICAL FUNCTIONALITY

Invacare® Rea® **Flex3**



Dual Stability System

Combined with the seat-tilt function, the **DSS** helps to lower or avoid friction and shear, which prevents Decubitus.

Rea® is the only company to offer **DSS**, which is fully Patented. **Rea®** was one of the first chairs designed to use innovative body length compensation (or virtual turning points) in backrests and legrests.

Weight shift in tilt. Maintains the stability of the chair as the person shifts their body weight forward. Additionally, the front part of the chair is raised by only 40-60 mm.

Body length compensation. The seat and the legrests adapt perfectly to the person's body length by simply changing the angle of the backrest and legrests.

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2012 Invacare International Sàrl.
All rights reserved.


All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Flex3 - UK - 03/2012 - 1553824



A model to suit you and your lifestyle

Suggested backrest options, plate or tension adjustable



Find the right wheelchair
 → *Rea® Azalea®* suitable for all people; available in attendant or occupant propelled variants



→ *Rea® Azalea® Tall* optimal comfort for taller people



Rea® Azalea® Max for heavier people



Rea® Azalea® Base for those who require more specialist seating systems



Rea® Azalea® Minor for shorter individuals or those of a lighter, slighter build

→ Best fits Flex3 system

Select the back rest system



Tension adjustable for people with specific seating requirements



Backrest plate for people requiring a particularly firm or stable backrest

Find the right support



Laguna aids lateral posture



Mistral provides targeted support in the waist area



Passad 2 provides support around shoulders and lower spine



The **lateral cover** combined with tension-adjustable backrest can be shaped to suit individual posture requirements



Shoulder High for tall people who require support in the regions of the shoulders and lower spine

All backrest cushions are available in plush or Dartex fabric.



Mix and Match	Backrest plate	Tension adjustable
Laguna	X	X
Passad 2	X	X
Shoulder high	X	X
Mistral 2	X	X
Lateral cover	n/a	X