

Invacare® | **Stream™**

NEW



Maintain independence

The **Stream**, a new power wheelchair from Invacare is a popular model suitable for those with reduced mobility. The **Invacare Stream** is ideal for both indoor and outdoor use. To ensure a comfortable sitting position and a relaxing ride, the seat can be adjusted to suit customers unique requirements. A dedicated lighting kit is available as an option, thus extending the possibilities of use for this model.



Yes, you can.®



Simple, straightforward and sturdy

The *Invacare Stream* offers a practical, functional design ideal for basic mobility. Thanks to 50 Ah batteries, this power wheelchair can achieve a range of 25 km. Simple settings for seat width, depth, height and angle are integrated into the wheelchair design and can be easily adjusted to suit particular requirements.

Ease of use

The *Stream* has been designed with ease of use in mind. Customers are free to travel with confidence and in comfort. The optional light kit is an added feature, providing enhanced safety for outdoor travel.

Features and Options



Light-kit

Available as an option - includes bright headlights with indicators.



Rear wheels

14" pneumatic rear wheels with studded profile.



Swing in/out footrest

The leg rests swing in and out for easy transfer and transportation.

Technical data

Stream	440 – 490 mm ¹	410, 460, 510 mm	425 – 500 mm ²	480 mm ²	290 - 360 mm ²	300 - 430 mm ²	up to 22°		
Stream	-10° - 30°	640 mm ³	1000 mm	1030 mm	105 kg				
Stream	130 kg	50 Ah	1660 mm	50 mm	11%				
Stream	2 x 180 W	25 km ⁴	6 or 9 km/h						

- 1/ seat cushion 43 cm
- 2/ measured without seat cushion / upholstery
- 3/ total width is depending on seat adjustment
- 4/ according to ISO 7176-4:2008

Invacare Limited
 Pencoed Technology Park - Pencoed -
 Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk
 Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County
 Dublin - Ireland
 Tel: +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.


Stream - UK - 10/2013

Frame colour



Gloss aluminium

Upholstery



Black nylon



More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.



Yes, you can.®