

# DOLOMITE | LEGACY



## The original

You want to have a reliable, good quality walker. A walker with more possibilities than a basic walker but without the most advanced adjustments and accessories.

**DOLOMITE LEGACY** is the original that offers you precisely this. That's why it has also been the inspiration for lots of similar walkers over the years.





**Walking area**

Large walking area inside the rollator frame is a signum for **DOLOMITE**. The flip-up seat gives premium walking space. Makes it easier to reach and access things. The rear forks gives a smooth inside, being friendly to toes and shoes.

Structured seat with anti-slip function, rounded rear edge and generous sitting area, means that the user can sit securely and comfortably for a long time.

**Safe when folding**

Folding handle positioned high-up so that the user avoids having to bend down and loose balance.

**Features and Options**



**Backrest**

Gives secure seating. Tray and basket helps transporting items.



**Kerb-climber**

Simplifies climbing kerbs or similar.



**Soft grips**

A softer alternative to anatomic grips. Find further options in separate leaflet.

**Technical data**

Legacy 520	60 cm	66 cm	46 cm	77 - 85 cm	54 cm
Legacy 600	60 cm	71 cm	46 cm	84 - 93 cm	62 cm

Legacy 520	40 cm	20 cm	7,1 kg	125 kg
Legacy 600	40 cm	20 cm	7,4 kg	125 kg

INVACARE® DOLOMITE AB  
 Växjövägen 303 75  
 Box 200, S-343 75 Diö  
 Tel. +46 (0)476 535 00  
 Fax +46 (0)476 535 99  
 www.dolomite.biz

© 2010 Invacare® International Sarl  
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Legacy - EU - 09/2010

**Colour**



RAL 5022

**Article code**

**Legacy 520**      12052  
**Legacy 600**      12050



**Yes, you can.®**