



Invacare H291 Cadiz

- Easy to handle seating height adjustment from 360 to 510 mm - in 13 steps
- Ergonomic seat surface
- Stable and secure stand due to angled rubber feet
- High load capacity: 143 kg



Invacare H296 Cadiz

- Easy to handle height adjustment from 360 to 510 mm - in 13 steps
- Ergonomic seat surface and ergonomic shaped backrest
- Stable and secure stand due to angled rubber feet
- Easy to clean


Technical data

	Invacare H291 Cadiz	Invacare H296 Cadiz
Shower stool w/o backrest	1544791	
Shower stool with backrest		1544829

Dimensions		
Seat width	510 mm	510 mm
Seat depth	310 mm	310 mm
Seat height	360 - 510 mm	360 - 510 mm
Backrest width	–	390 mm
Backrest height	–	230 mm
Total width	440 mm	440 mm
Total depth	415 - 460 mm	415 - 460 mm
Total height	405 - 555 mm	790 - 940 mm

More info		
Load capacity	143 kg	143 kg
Weight	3 kg	3 kg
Colour	White	White

Invacare H291 Cadiz



Invacare H296 Cadiz

