

Invacare® | **Propad Leg Trough**



Shaped leg support - provides comfort and promotes recovery

The *Invacare Propad Leg Trough* has been designed in consultation with clinicians to provide maximum support for the whole leg. The **Leg Trough** elevates the heels and allows them to be suspended over the edge of the trough. This provides maximum pressure relief for the vulnerable heel area.

Propad



Yes, you can.®



Features and Options

Provides comfort and support to promote recovery

The leg trough offers optimum support and protection for the whole leg, reducing pressure and promoting recovery.

Relieves pressure on vulnerable heels

Elevates the heels and allows them to be suspended over the edge of the trough, providing pressure relief for the vulnerable heel area.



Heavy duty water-resistant cover

A two-way stretch, vapour permeable and anti-static polyurethane (PU) water-resistant cover.



Cover includes a zip and is removable

The cover can be unzipped and removed from the foam core to allow laundering of the cover.



Adjustable side-squeeze clip buckles

The adjustable buckles enable appropriate and comfortable positioning of the leg.

Technical data

Propad Leg Trough



Length:
Small: 510mm
Medium: 610mm
Large: 660mm



Width:
110 mm



Height:
115 mm



Weight:
0.8 kg



Cleaning: 80°

Manufactured to European Standards:
Interior Foam: BS 5852: 1982 Crib 5, BSEN597-1,2
Cover: BS7175: 1989 Crib 5

Colour



Blue

Warranty 1 Year

