

# Invacare® | **Dragon**®



## **Compact and Versatile**

The *Invacare Dragon* is a solid and cost effective power chair that provides excellent indoor and outdoor manoeuvrability. The new seat height and tilt feature can easily be adapted to suit the needs of the client and can be achieved without taking the seat apart.



**Yes, you can.®**



**Features and Options**

**Mobility and Comfort**

The robust frame is available with a choice of seating that makes the *Invacare Dragon* a reliable wheelchair for daily use. The compact dimensions provide good manoeuvrability and the spring suspension of the motors absorbs shocks for increased comfort and the large 60 ah batteries provide a longer range. The powerful **Shark II** controller enhances driving with its fingertip precision and is programmable to suit individual needs.

**Adjustability and Easy Service**

*Invacare Dragon* features a range of seating, a manual settable seat tilt and 6 km/h motors. It is adjustable to suit the needs of clients and all parts of the power chair are accessible for ease of servicing.



**Tilt and Seat Height**

With this new feature seat height and tilt can be easily adapted to the clients needs.



**Contour Seat**

Contoured upholstered seat for improved comfort.



**Shark II Control**

The version of the **Shark** remote offers the possibility to control lights with indicators and up to two electric adjustments.

**Technical data**



Standard Seat	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	480 / 540 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>
Contour Seat	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	560 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>
Comfort Seat	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	560 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>



Standard Seat	-10° - 30° <sup>3</sup>	610 mm <sup>4</sup>	980 / 1040 mm	1030 mm <sup>5</sup>	850 mm	100 kg <sup>8</sup>	140 kg
Contour Seat	-10° - 30° <sup>3</sup>	610 mm <sup>4</sup>	1060 mm	1030 mm <sup>5</sup>	850 mm	105 kg <sup>8</sup>	140 kg
Comfort Seat	-10° - 30° <sup>3</sup>	610 mm <sup>4</sup>	1060 mm	1030 mm <sup>5</sup>	850 mm	105 kg <sup>8</sup>	140 kg



2 x 12 V 60 Ah	1660 mm	60 mm <sup>6</sup>	15% (8.5°)	2 x 180 W	30 km <sup>7</sup>	6 km/h
----------------	---------	--------------------	------------	-----------	--------------------	--------

INVACARE Limited  
 Pencoed Technology Park - Pencoed -  
 Bridgend CF35 5AQ  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)  
 Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)

INVACARE Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County  
 Dublin  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2011 Invacare International Sàrl  
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Dragon - UK - 07/2011 - 1539041

**Frame colour**



Blue metallic

- 1/ seat cushion 430 mm
- 2/ varies with seat height
- 3/ in 6 steps, 5° -30° with gas spring, option
- 4/ depending on seat adjustment
- 5/ incl. leg rest and antitippers
- 6/ 100 mm with kerb climber
- 7/ acc. to ISO 7176-4:2008
- 8/ total weight is depending on configuration

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.



**Yes, you can.<sup>®</sup>**