

Invacare® | *Dragon® Vertic*



Access all areas – with the push of a button

The *Invacare Dragon Vertic* is an agile power chair that not only provides mobility to an active customer, but will bring you face to face with other people. Its comfortable seating system easily moves you into a standing position by a simple joystick command, providing increased access and excellent mobility.



Yes, you can.®



Mobility when upright

The *Invacare Dragon Vertic* greatly enhances activity in your daily life with the vertical function. With mini front castors for extra stability you can even drive whilst in the upright position giving you increased access and mobility. The process of going upright is always from the seating position (driving) into vertical. The speed for moving up and down can be controlled via the joystick.

Easy transfer

The knee fixation can be swivelled left or right, or be taken away entirely for easier transfer. The transfer is also facilitated by swivelling-up armrests and foot plates.

Comfortable positioning

Seat depth, legrest length, backrest tension and height can be adapted to the needs of the customer, and the angle of

the backrest can be electrically adjusted for optimum comfort when driving and when upright.

Features and Options



The armrest can be turned upward for an easier transfer.

The headrest is individually adjustable in width, height and horizontal position.

The sliding backrest minimises the shearing forces and is height adjustable.

Features and Options



The hip support is width adjustable and removable for transfer.



The adjustable back rest angle optimises comfort in seating and upright position.



The foot plates can be swivelled up for easy transfer.



Puncture proof rear tyres available.



Mini castors allow driving in the upright position.



Tray
Convenient tray for work and for trunk support, can be swivelled down.



Knee support
The knee support is adjustable in depth, width and height it can be swivelled left and right, or can be taken off entirely



Shark II
Easy steering with **Shark II**. Variable speed for vertical function and for backrest angle adjustment.



On the road version
Front and rear lights provide safety when driving on public roads. The light is activated via the Shark controller.



Technical data



Dragon Vertic 38,5 - 46 cm
45 cm incl. armpads



38 - 46 cm



50 / 55 cm



60 / 65 cm



26 - 35 cm



31 - 45 cm



167 / 172 cm



Dragon Vertic 6° - 22°



48.5 cm¹



116 / 121 cm



105 cm²



85 cm



131 kg³



100 kg



Dragon Vertic 2 x 12V 60 Ah



166 cm



4 cm



17,6 %



2 x 180 W



25 km⁴



6 km/h

- 1/ total width is depending on seat adjustment
2/ incl. legrest and anti-tippers
3/ depending on configuration
4/ acc. to ISO 7176-4:2008

Frame colours



Blue metallic

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© Invacare® International Sarl 2010
All rights reserved.

All measurement have been done under the Invacare® standard test conditions, unless quoted differently. All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Dragon Vertic - EU - 08/2010



Yes, you can.®

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.