



PRODUCT INFORMATION UPDATE - January 2008

Power Wheelchairs - Mirage

Bulletin PU240
Date: 16/01/2008

Please find below details of the latest product changes for our Powerchair range.

Item	Description of change	New Part No. (if applicable)	Introduction date	Picture
MIRAGE	In order to reduce the amount of positive camber on the rear wheels of the Mirage, we have removed the two motor spacing strips (part number 1441980-0216)		Jan-08	
	<div style="border: 1px solid black; padding: 5px;">The serial number of first chair with this change is: 060812178140</div>			

If you require any additional information regarding any of the above changes, please contact your local

CUSTOMER SERVICES

TEL: 01656 647327

FAX: 01656 649016

EMAIL: uk@invacare.com