

Invacare® | **Harrier XHD**



**Power and strength**

The *Invacare Harrier XHD* is a sturdy well designed heavy duty chair offering reliability and durability for the active person. With all the features you expect from the *Harrier XHD* offers comfort and support for day long use.



**Yes, you can.®**



**Features and Options**

**Comfort and Security**

The *Invacare Harrier XHD* is a classic powered wheelchair with many features as standard such as firm upholstery, adjustable armrests and puncture proof tyres. With two 50 Ah batteries you can be assured that the *Harrier XHD* has the extra power and strength to cope efficiently outdoors. A range of control options is available to support the individual's needs.

**Travel in style!**

The *Harrier XHD* features the padded seating system with firm back and deep padded seat providing a fully supportive and comfortable sitting environment for day long use. With a range of options and accessories the *Harrier XHD* provides a very flexible power chair with a carrying capacity of up to 160 kg. Simply plug into any household mains supply to recharge the batteries.



**Firm upholstery**  
A fully padded seat and backrest offer additional support and day-long comfort.



**Padded side guards**  
Offer extra comfort and support.



**Swing away footrests**  
For easier access to the front of the chair when getting in and out.

**Technical data**

<b>Harrier XHD</b>	460 - 510 - 560 mm	370 - 500 mm	580 mm	420 - 570 mm	180 - 260 mm <sup>1</sup>	420 - 510 mm
<b>Harrier XHD</b>	0° - +70°	24°	660 - 760 mm	960 mm	typ. 1130 mm	820 mm
<b>Harrier XHD</b>	76 kg <sup>2</sup>	max. 160 kg	50 Ah	1940 mm	max. 100 mm	9°
			Given values may vary with configurations			
<b>Harrier XHD</b>	75 W	20 km <sup>3</sup>	<sup>1/</sup> according to ISO 7176-4:2008 <sup>2/</sup> total weight is depending on configuration <sup>3/</sup> The maximum range is dependent upon several factors, which include environment and terrain, and battery conditions etc.			

INVACARE Limited  
Pencoed Technology Park - Pencoed -  
Bridgend CF35 5AQ  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: uk@invacare.com  
www.invacare.co.uk  
Sales Order E-mail:  
ordersuk@invacare.com

INVACARE Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: Ireland@invacare.com  
www.invacare.ie

© 2011 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare reserves the right to alter product specifications without prior consultation.

Harrier XHD - UK - 07/2011  
1441996

**Frame colour**



Black

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.



**Yes, you can.®**

