

# Invacare® | **Harrier Plus**



## **Practical and reliable**

The *Invacare Harrier Plus* is a classic powered wheelchair with a range of features that give the ability to get about in comfort and safety. This flexible chair is suitable for both adults and children, and is fully adjustable. Enjoy life to the full, indoors and out with the *Harrier Plus* - the ideal choice for the more active person who looks for reliability at all times.



**Yes, you can.®**



### Flexibility, choice and safety

The *Invacare Harrier Plus* has a standard canvas seat and back, both of which are slightly padded for comfort, offering gentle support combined with practicality. The *Harrier Plus* is also available with an optional firm back and deep padded seat, giving a fully supportive, comfortable sitting environment for day long use. The *Harrier Plus* has a range of features for added safety and comfort including parking brakes for use when stationary, a standard transportation kit, and a central kerb climber for navigating the ups and downs of the high street!

### Choice and versatility

Available in a range of seat sizes, the *Harrier Plus* is ideal for growing youngsters, with a choice of backposts and a range of adjustable footrests which allow the chair to be easily and economically upgraded. This flexibility means that it is possible to 'tailor-make' the *Harrier Plus* to suit individual needs.

### Added value in a classic chair

Additional standard features such as the puncture proof tyres and castors combined with top quality construction make the *Harrier Plus* the ideal choice in terms of durability and reliability. The maintenance-free batteries give a range of up to 25 km on one charge.

The *Harrier Plus* is the smart and practical choice for the more active person - young or young-at-heart - in today's busy world!

## Features and Options



#### Folding backrest

Reduces the overall size of the chair when folded for transportation or storage.



#### Chin control

For those with very limited movement.



#### Electrically powered reclining backrest

Gently reclines for a relaxed position making it possible to rest in the chair.





## Features and Options



**Upholstered side and arm pads**  
Giving extra comfort and support.



**Elevating legrests**  
A choice of manual or electric elevating legrests provide support and relief when resting.



**Footrests**  
Fully adjustable footrests gives the correct leg angle.



**Firm headrest and extended height backrest**  
For added support and comfort.



**Backrests**  
A choice of backrest options to reduce or increase the seat depth as the person grows.

## Technical data

<b>Harrier Plus</b>	330 - 460 mm <sup>1</sup>	300 - 500 mm	500 mm	410 mm	210 - 350 mm	80 - 190 mm	0° - +114°
<b>Harrier Plus</b>	4° - 12°	550 - 625 mm	920 mm	1060 mm	720 mm	65 kg <sup>2</sup>	max. 114 kg
<b>Harrier Plus</b>	50 Ah <sup>1</sup>	1920 mm	max. 100 mm	10°	170 W	25 km <sup>3</sup>	max. 6,5 km/hr

- 1/ 50Ah batteries as standard except for seat width 380mm which has 30Ah batteries  
 2/ according to configuration  
 3/ The maximum range is dependent upon several factors, which include environment and terrain, and battery conditions etc.  
 Given values may vary with configurations

Range: These results are achieved under test conditions stipulated to ISO 7176-4:2008.

## Frame colour



Black

Please notice colors may vary slightly from those displayed above.

INVACARE Limited  
 Pencoed Technology Park - Pencoed -  
 Bridgend CF35 5AQ  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: uk@invacare.com  
 www.invacare.co.uk  
 Sales Order E-mail:  
 ordersuk@invacare.com

INVACARE Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County  
 Dublin  
 Tel: +353 1 8107084  
 Fax: +353 1 8107085  
 Email: Ireland@invacare.com  
 www.invacare.ie

© 2011 Invacare International Sàrl  
 All rights reserved.

All information quoted is believed to be correct at time of print.  
 Invacare reserves the right to alter product specifications without prior consultation.

Harrier Plus - UK - 07/2011 - Q99621