

Invacare® | *Spectra Blitz*



Go and grow with the *Invacare Spectra Blitz*

The *Invacare Spectra Blitz* is a versatile power chair for children. Its compact design provides excellent mobility and driving performance. It features a highly adaptable seat and allows precise adjustment for optimal comfort and posture management.



Yes, you can.®



Features and Options



Height adjustable back rest.



Adjustable seat width from 250mm to 390mm.




Rear wheel suspension for better shock absorption.

Technical data

| | | | | | | | |
|----------------------|--------------|--------------------|--------|--|-----------|--------|--------|
| | | | | | | | |
| Spectra Blitz | 250 - 390 mm | 250 - 380 mm | 450 mm | 370 - 440 mm | -5° - 25° | 560 mm | 940 mm |
| | | | | | | | |
| Spectra Blitz | 910 mm | 1240 mm | 50 mm | 14.1% (8°) ³ | 2 x 40 Ah | | |
| | | | | | | | |
| Spectra Blitz | 2 x 106 W | 21 km ¹ | 6 km/h | 79,5 kg (incl. batteries) ² | 75 kg | | |

- 1/ according to ISO 7176-4:2008
- 2/ depending on configuration
- 3/ The maximum safe slope depends on the configuration. For further details please see the user manual.

Frame colours



INVACARE® Limited
 Pencoed Technology Park - Pencoed -
 Bridgend CF35 5HZ
 Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
www.invacare.co.uk
 Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County
 Dublin
 Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare® International Sarl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare® reserves the right to alter product specifications without prior consultation.

Spectra Blitz - UK - 02/2011 - 1450089



Yes, you can.®



More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.