

Invacare® | **Vicair® Academy Junior Range**



**Warranty
2 Years**



The optimal solution in pressure reduction and positioning for children

The *Invacare Vicair Junior range* provides a unique combination of pressure re-distribution, comfort and optimal positioning for children, making it suitable for those at 'Very High Risk'* of developing pressure ulcers. The comfortable back and seat cushion can be adjusted in cell density and customised to fit the child's body contours, offering high quality support. The **SoftCells™** provide a dynamic seating surface that can be adjusted to accommodate a range of postural asymmetries.



Yes, you can.®



The Vicair Junior Range

Unique design

The clinically proven cushions offer a unique and innovative approach to postural support for children. The range contains small air filled, pyramid shaped **SoftCells™** in strategically placed compartments. These flexible, fabric coated cells have a low friction level when in contact with each other, helping to reduce shear levels.

Light weight and low maintenance

The range is made up of a paediatric cushion and back system and aims to provide a high level of comfort and stability for children whilst being safe, easy to maintain and lightweight. There is no risk of deflation. If a **SoftCell™** does get damaged, the individual cell can be replaced. There is little staff intervention required as daily checks are not necessary plus it is ultra lightweight making it easy to lift and transfer.



Academy Junior Cushion





The **Academy Junior** cushion is based on the **Academy Vector**** but has been specifically designed for children. With six compartments, the cushion provides optimal pressure re-distribution and stability and offers a comfortable surface for children to sit on. The cushion is ideal for children who need additional pelvic stability.

Academy Junior Back

The Back system consist of a cushion with 8 compartments and a removable ABS plastic back shell, providing comfort and good lateral support. Supplied with a velcro attachment, the system can be fitted to the back posts of most types of paediatric wheelchairs.



Technical Data

				
Seat cushion	Size (w x d) 250x250mm 300x300mm 350x350mm	Weight: 0.4kg	Max User Weight: 15.5 stone/100kg	Outer cover: washable to 60°. Inner cover: Hand wash. SmartCells: Hand wash
Back	Size: 250, 300, 350mm	0.65kg	11 stone/70kg	

Cover options

Comfair: Multi-stretch breathable top fabric allows air to circulate, provides comfort and reduces the risk of moisture/ heat build up.

Incontinence: Multi-stretch water resistant and vapour permeable fabric with an easy to maintain wipedown surface.

* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

** For information on the **Academy Vector** cushion, please refer to the *Invacare Vicair Academy* range brochure.

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2012 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Vicair Junior - UK - 03/2012 - 1554355



Yes, you can.®

