

Part	Feature	Benefit
Frame	3 different rear bracket positions. Standard, active, amputee 15 rearwheel positions	Adjustable according to user need and activity level
	Maximum user weight: 125 kg	A reliable product meeting the need of the majority of the population
	Overall width is seat width + 19 cm	Enables the user to pass through even small doorways
	Rigidity	Frame construction
	Driving Performance	Increased by 15 rearwheel positions Increased by construction of rear bracket. Low weight 14,5 kg
Seat	Seat angle can be adjusted in 3 angles 0, 3 and 6 degree	Allows the user to obtain correct seating position
	Seat height from 38,5 cm to 51 cm	Allows even a very small user or a hemiplegic person an individual assessment
Backrest	Standard 7 degree or 0 degree Angle adjustable -6 / +12 degree Reclining Backrest 0 – 30 degree Comfort backrest	Insures correct seating for the individual need of the end user
Armrest	24 cm standard height of armrest long and short. Increased length behind the backrest tube. Height adjustable from 18 cm – 28 cm Comfort version with depth adjustment	Insures correct position highly improved when using a reclining backrest because the armpad is longer.
	Armrest adjustable in height	Allows the user to obtain the correct support while seating and driving
	Flip up armrest	Easy side transfer in and out of the wheelchair
Footrest	Release mechanism manoeuvres 3 ways Left, right and pulling to the front Very stable footrest Interchangeable left and right side Good position of the calf straps or heel straps	Insures the user and assistant will experience high security and correct seating.

<p>Footplates</p>	<p>Fixed Footplates Adjustable in depth 5 cm Angle offer infinity adjustment Adjustable in steps of 1 cm Adjustable in length form 32 – 44 cm and 41 – 53 cm with 2 different inner tubes Step adjustment insures footrest in position</p>	<p>Footplates and adjustments insures that the end user will be seated correctly in the wheelchair with good positioning.</p>
-------------------	--	---