

Invacare® | **Phoenix® Plus**



A great deal more for less

Invacare has taken the original *Invacare Phoenix* and re-styled and updated it to produce a more stylish, more flexible and comfortable power chair, while retaining the core features of easy to use controls, adjustability and quick and easy dismantling. Last but not least, the *Invacare Phoenix Plus* is now available at an even more competitive price.



Yes, you can.®



The smarter choice

A practical, compact power wheelchair, the **Invacare Phoenix Plus** is suitable for both indoors and outdoors. The **Phoenix Plus** comes with 40 Ah batteries as standard, increasing your mobility and independence. The additional features such as angle adjustable footplates, padded clothes guards and armrests improve all-round comfort and positioning, minimising fatigue and tiredness.

Practicality

The **Phoenix Plus** is the ideal choice for manoeuvrability and versatility with its compact footprint.

This practical wheelchair can be easily dismantled for transportation or storage.

Features and Options



Folding Back

Reduces overall size when folded for transport or storage.



Height Adjustable Armrest

The armrests are easy to remove and adjustable in height.



Alloy wheels

Durable alloy rims to ensure a trouble-free ride.

Technical data



Phoenix Plus 400 / 440 / 480 mm



400 - 440 mm



500 mm



440 (700) mm



220 - 320 mm



360 - 460 mm



2°

0° (30°)



Phoenix Plus

570 mm¹



970 / 1230 mm



1060 mm



730 mm



68 kg²



114 kg



Phoenix Plus

2 x 12V / 40 Ah



1940 mm



50 mm



max 8.5°



2 x 110 W



20 km³

1/ depending on seat adjustment

2/ total weight is depending on configuration

3/ Please note that range quoted is based upon ISO 7176-4:2008 and may vary dependent upon many factors (user weight, terrain, battery condition etc)

Please note colours may vary slightly from those displayed above.

Frame colours



Black



Transparent Purple

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.

INVACARE® Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Phoenix Plus - UK - 09/2010
Q99642



Yes, you can.®