

Invacare Ireland Ltd  
Unit 5, Seatown Business Campus,  
Seatown Road,  
Swords, County Dublin  
Tel: (00353) 1 8107084 Fax: (00353) 1 8107085  
email: [ordersireland@invacare.com](mailto:ordersireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)



**Yes, you can:**

Powered Wheelchairs  
Ireland Prescription Form  
LPF1U2PHOENIXPLUS150513IRE

**July 2013**  
**IREPHNX3W**

Dealer Name:	Account No :
Address :	Customer Name :
Tel No:	Order No :
Delivery address :	Order Qty :
	Order Date :

## **Invacare® Phoenix Plus 3**

**INVACARE IRELAND LTD**

**Max User Weight 114kg (18 Stone)**

**Product P/N :- 1552902** ☐



### **Standard Features**

Black Nylon Upholstery  
Sling Seat  
45cm Fixed Backrest Height  
48cm Seat Width  
44cm Seat Depth  
Swingaway Detachable Legrests  
Fixed Footplates  
Moulded Armrests and Clothesguard - Wide  
Penny and Giles VSI Controls R/H  
Swingaway J/Stick Bracket R/H  
Black Frame Colour  
200 x 50mm Microcellular Castors  
315mm Alloy Split Rim Solid Infil Rear Wheels  
40 amp Batteries  
Charger  
Central Kerb Climber

**Base Price €1,917**