

DOLOMITE | *ALPHA*



For users who require extra support

DOLOMITE ALPHA is a family of solutions giving added support for gait training, arthritis, orthopaedic needs, after surgery and ...

ALPHA does not only offer an excellent support, when other products fail to give the right arm positioning or functioning handles **DOLOMITE ALPHA** is the product to try.



Infinite adjustable handle

The basic version of **DOLOMITE ALPHA** has all the necessary height, width and depth adjustments that are needed to provide the user who requires extra support with a good walking position. In addition it has rounded armrests that give comfortable support for the forearms.

Ready for advanced customisation

DOLOMITE ALPHA ADVANCED - is unique in the opportunities it presents for multi-directional adjustments thanks to its ball- and-socket design. This enables the Advanced model to live up to the extremely high demands that patients who are in pain place on handle adjustments. The angle of the handles can be set as required, as can the brakes operating direction and height of the brake handle.

Features and Options



Platform support

Large rest area with platform support.



Adapt handle with one-hand-brake

For stroke users.



One-hand-brake

Operate both brakes with One hand.

Technical data



Basic	65 cm	70 cm	51 cm	95 - 124 cm	62 cm
Advanced	65 cm	70 cm	51 cm	95 - 124 cm	62 cm
Basic Platform	65 cm	70 cm	51 cm	103 - 124 cm	62 cm
Advanced Platform	65 cm	70 cm	51 cm	103 - 124 cm	62 cm



Basic	43 cm	20 cm	10,9 kg	125 kg
Advanced	43 cm	20 cm	11,4 kg	125 kg
Basic Platform	9 - 27 cm	20 cm	12,4 kg	125 kg
Advanced Platform	10 - 30 cm	20 cm	12,9 kg	125 kg

Colour



RAL 5022

Art. nr.

Basic	12230
Advanced	12220
Basic Platform	12230-02
Advanced Platform	12220-02

INVACARE® DOLOMITE AB
Växjövägen 303 75
Box 200, S-343 75 Diö
Tel. +46 (0)476 535 00
Fax +46 (0)476 535 99
www.dolomite.biz

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Dolomite Alpha - EU - 09/2010



Yes, you can.®