

DOLOMITE | *MAXI+*



Stronger and smarter

DOLOMITE MAXI+ is developed for groups of users needing a high weight capacity and a generous frame size. The new **DOLOMITE MAXI+** has a reinforced frame with a maximum weight capacity of 200 kg/440 lbs.





Extra space

The wide seat and the additional 10 cm space between the handles allows for comfortable seating.

The seat can be folded up to gain more space for walking and reaching things, it also facilitates to opening doors.

Adaptable to individual needs

DOLOMITE MAXI+ is available in two different heights. The long height adjustment has an overlapping adjustment range. Choose the higher seat height when having difficulties to change from sitting to standing.

DOLOMITE MAXI+ has a wide range of accessories to adapt the product further to individual needs.

Features and Options



Comfortable backrest



Slow-down brake

A safety feature preventing the walker from rolling away.



Cane holder

Makes the cane easily accessible without getting in the way.

Technical data

Maxi 55	69 cm	65 cm	53 cm	72 - 89 cm	56 cm
Maxi 65	69 cm	74 cm	53 cm	81 - 98 cm	65 cm

Maxi 55	47 cm	20 cm	9,6 kg	200 kg	
Maxi 65	47 cm	20 cm	9,9 kg	200 kg	

Colour



RAL 5023

Article code

Maxi 55 12122
Maxi 65 12120

INVACARE® DOLOMITE AB
Växjövägen 303 75
Box 200, S-343 75 Diö
Tel. +46 (0)476 535 00
Fax +46 (0)476 535 99
www.dolomite.biz

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Maxi+ - EU - 09/2010



Yes, you can.®

