

Invacare® | **Ben 9 Plus**



Simple and practical

The *Invacare Ben 9 Plus* is the entry level chair in the Invacare manual range. The lightweight construction of this transit chair makes it easy for even an elderly carer to propel providing a comfortable and adaptable chair.





Features and Options

Added value

Invacare has built in many beneficial features as standard on the **Invacare Ben 9 Plus**, all designed to make your life easier. These include the hinged backrest and detachable armrests for ease of transfer and to reduce the overall size when folded for transportation. In addition, the footrests are swing away and detachable to aid the transfer process.

Smart and practical

The **Ben 9 Plus** is available with a black frame, with a Black nylon canvas seat and backrest.

A wide range of options are available, from contoured cushions to extended backrests and headrests for added support and comfort. Simply choose the options to suit your individual needs.



Folding backrest

Reduces overall size when folded for transportation or storage.



Swing away footrests

Footplates fold up as shown, and footrests swing away, allowing easier access in or out of the chair.



Foldable cross frame

Lightweight frame folds flat for easier portability.

Technical data

38 - 45 cm	40,6 - 43,2 cm	48 cm	43 cm	22 cm	33 - 43 cm

5°	8°	Seat width + 13 cm	96,5 cm	91 cm	67 cm

17 kg	127 kg	Crashtested ISO 7176-19*

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Ben 9 Plus - EU - 08/2010

Frame colour



Black

Please note colours may vary slightly from those displayed above

