



Yes, you can:



e-spare cat Now Available Online
 at www.invacare.ie

Customer Name :	Account No :
Address :	Order No :
Tel Number :	Order Qty :
Delivery Address :	Order Date :

Invacare® Ben 9NG HD Stock Order Form
 INVACARE IRELAND LTD
 Max User Weight 160Kg (25 stone)

- Standard Features:** CE
- Transit (Attendant) Propulsion or Self (Occupant) Propulsion
 - Seat Width Range from 510mm to 610mm (20" to 24")
 - Seat Depth either 430mm (17") or 450mm (18")
 - 485mm Seat Height (19")
 - Half Folding 10° Backrest, 420mm Backrest Height (16.5")
 - Black Nylon Upholstery
 - Full Length, Fixed Height, Removable Armrest
 - Swing Away and Detachable Legrest Hangers
 - Length Adjustable, Flip Up, Plastic Footplates
 - Solid Castor and Rear Wheel Tyres
 - Composite 315mm (12") Rear Wheels MCP Tyres
 - or 610mm (24") Self Propel (Spoked, Solid Tyres)
 - Lapbelt (adjustable both sides)
 - Push To Lock Brakes

Invacare® Ben 9NG HD Transit

(Attendant) Propulsion c/w Composite 315mm (12") Rear Wheels MCP Tyres

Bulk order quantity (Please state relevant order quantity next to the product configuration required). *NB: All sizes are approximate*

Standard Net Price €305

1550174	510mm x 430mm (20" x 17"), Black, Transit	N/A	_____
1550175	510mm x 460mm (20" x 18"), Black, Transit	€33	_____
1550176	560mm x 430mm (22" x 17"), Black, Transit	€33	_____
1550177	560mm x 460mm (22" x 18"), Black, Transit	€65	_____
1550178	610mm x 430mm (24" x 17"), Black, Transit	€33	_____
1550179	610mm x 460mm (24" x 18"), Black, Transit	€65	_____

Invacare® Ben 9NG HD Self Propel

(Occupant) Propulsion c/w 610mm (24") Self Propel (Spoked, Solid Tyres)

Bulk order quantity (Please state relevant order quantity next to the product configuration required). *NB: All sizes are approximate*

Standard Net Price €305

1550184	510mm x 430mm (20" x 17"), Black, Self Propel	N/A	_____
1550185	510mm x 460mm (20" x 18"), Black, Self Propel	€33	_____
1550186	560mm x 430mm (22" x 17"), Black, Self Propel	€33	_____
1550187	560mm x 460mm (22" x 18"), Black, Self Propel	€65	_____
1550188	610mm x 430mm (24" x 17"), Black, Self Propel	€33	_____
1550189	610mm x 460mm (24" x 18"), Black, Self Propel	€65	_____