

Invacare® | **Duo / Duo HD**



Versatile and affordable

The *Invacare Duo* offers an economical solution to your mobility needs. With its wide range of seat widths 15" - 18" maximum occupant weight 20 stone and 20" - 24" maximum occupant weight 25 stone (heavy duty version) which makes for a comfortable and adaptable chair.



Yes, you can.®



Adaptable

This practical wheelchair has been developed from the **Ben 8L** and offers many features as standard to meet client requirements. It's steel construction with 22" quick release wheels and folding backrest makes for a compact and manoeuvrable chair.

The *Invacare Duo* features a rear wheel set back option for increased stability; this can be easily achieved by rotating the rear axle block.

Practical

The chair can be easily adapted to accommodate alternative seating options, depending upon client requirements. A simple cross frame design with swing away footrests and removable armrests allows for easier access in or out of the chair. For ease of storage and transportation the chair will dismantle and fold away.

Features and Options



Rear wheel forward position



Removable armrests



Easier side transfers

Technical data

Invacare Duo / Duo HD	40 - 62,5 cm	42 - 47 cm	48 cm	47 cm	22 cm
					35 - 44 cm

Invacare Duo / Duo HD	3,5°	10°	Seat width + 20 cm	68 - 72 cm (Folded)	104 cm
					78 cm

Invacare Duo / Duo HD	19,3 - 19,7 kg	125 / 160 kg (HD Version)
		Crashtested ISO 7176-19*

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Duo / Duo HD - EU - 08/2010

Frame colour



Black

