



Yes, you can.®

Invacare® Essential Basic Cushion


High quality at a cost effective price

The *Invacare® Essential Basic Cushion* is manufactured from highly resilient foam, providing a good level of comfort at a cost effective price, making this product ideal as an entry level wheelchair cushion.



Cushion Depths

The **Essential** is available in two depths, a 2" depth known as the **Essential Low Profile** and a 3" known as the **Essential Basic**.



Features




Ingress resistant cover as standard
Black smooth, high quality, ingress repellent zipped cover as standard. The fabric allows water to 'roll off', protecting clients clothing and the foam core from soiling. A plaiton interliner is also available as an option.



Available in a wide range of sizes
To view the cushion size range, please visit www.invacare.co.uk.


Technical data




Height:
Low profile 2" / 5.08 cm
Basic 3" / 7.62 cm



Max Client Weight:
3" 17 stone (108 kg)
2" 12 stone (76 kg)



Weight: 0.8 kg
Weight: 1 kg



Wide range of sizes



Cleaning: 80°C
Refer to instructions for use and relevant infection control guidelines.

Cover:
Ingress Resistant
Vapour Permeable

Cover Option



Black Smooth
(Standard)

Manufactured to European Standards:
BS7176-16



Yes, you can.®

Invacare® Essential Visco Cushion


Dual layered comfort

The *Invacare® Essential Visco Cushion* has a dual layer, using a combination of high density temperature sensitive visco elastic foam in conjunction with a highly resilient foam base to provide added support and longevity. The *Essential Visco Cushion* is suitable for those who require added comfort when in a wheelchair.



➤ **Visco elastic (memory) foam top**

The top layer of the cushion is made from high density visco elastic (memory) foam, which conforms to the shape of the body combining exceptional comfort with a fully supported feeling.



➤ **Robust Foam Base**

The high density foam base provides added support whilst preventing the cushion from bottoming out, ensuring longevity of the product.

Features



Ingress resistant cover as standard

Black smooth, high quality, water repellent zipped cover as standard. The fabric allows water to 'roll off', protecting clients clothing and the foam core from soiling. A plaiton interliner is also available as an option.




Available in a wide range of sizes


To view the cushion size range, please visit www.invacare.co.uk.

Technical data

Warranty
2 Years



Wide range of sizes



Height:
7.62 cm



Max Client Weight:
17 stone (108kg)



Cleaning: 80°C
Refer to instructions for use and relevant infection control guidelines.



Weight: 1.4 kg

Cover:
Ingress Resistant
Vapour Permeable

Cover Option



Black Smooth
(Standard)

Manufactured to European Standards:
BS7176-16