

# Aquatec® | **Sorrento Shower Collection**

NEW



## **The innovative, modern and modular choice for your bathroom**

The *Aquatec Sorrento Shower Chair and Stools* are part of the outstanding modular and modern products of the *Aquatec Sorrento range*. This range will convince you not only with great function but also with great style and design.

The non-slipping soft touch surfaces as well as the soft backrest provide a secure and very comfortable seat for the user. The height adjustable, folding legs as well as the fast and easy assembly cover various users' needs.

The all-corrosion free, modular and knock-down shower chair range is the best choice for users who need assistance in the bathroom for their personal hygiene.

**AQUATEC®**



**Yes, you can.®**



### Non-slipping Soft Touch Surface

The Aquatec **Sorrento** is the first shower chair with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.

### Height adjustable and folding legs

The Aquatec **Sorrento range** can be tool less adjusted to individual body sizes. With its folding mechanism it is ideal to take on vacation or for easy storage. Backrest and armrests are fast and easy to take off.

## Features and Options



Sorrento as Shower Chair



Sorrento as Shower Stool with armrests



Sorrento as Shower Stool



All versions also available with closed seat



Sorrento – looks great in every bathroom

## Technical data

	Sorrento Shower Chair	Sorrento Shower Stool w/ armrests	Sorrento Shower Stool	Sorrento Shower Chair w/o recess	Sorrento Shower Stool w/ armrests, w/o recess	Sorrento Shower Stool, w/o recess
Shower Chair	1532633			1532730		
Shower Stool w/armrests		1532660			1532729	
Shower Stool			1532661			1532728

Measurements	Sorrento Shower Chair	Sorrento Shower Stool w/ armrests	Sorrento Shower Stool	Sorrento Shower Chair w/o recess	Sorrento Shower Stool w/ armrests, w/o recess	Sorrento Shower Stool, w/o recess
Seat width (mm)	400	400	400	400	400	400
Seat depth (mm)	400	400	400	400	400	400
Seat height (mm)	460 - 610	460 - 610	460 - 610	460 - 610	460 - 610	460 - 610
Total width (mm)	520	520	520	520	520	520
Total depth (mm)	430	430	430	430	430	430
Total height (mm)	850 - 990	850 - 990	850 - 990	850 - 990	850 - 990	850 - 990
Range seating height adjustment (mm)	150	150	150	150	150	150
Width between armrests (mm)	450	450		450	450	
Armrest length (mm)	280	280		280	280	
Armrest height (mm)	195	195		195	195	
Backrest width (mm)	500			500		
Backrest height (mm)	380			380		
Backrest cover height (mm)	250			250		
Width hygiene recess (mm)	140	140	140			
Depth hygiene recess (mm)	170	170	170			
High load capacity	135 kg	135 kg	135 kg	135 kg	135 kg	135 kg
Weight	5,5 kg	4 kg	3,5 kg	5,5 kg	4 kg	3,5 kg
Colour						

Invacare International Sarl  
Route de Cité-Ouest 2 - 1196 Gland  
Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
hqeurope@invacare.com  
www.invacare.eu.com

© 2011 Invacare International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare® reserves the right to alter product specifications without prior consultation.

Sorrento Shower Collection - EU - 03/2011



Yes, you can.®

