

**AQUATEC®**

# Ocean E-Vip

**One carer operation\***



Taking care of the carer



With its full powered adjustment, the *Aquatec Ocean E-VIP* shower chair commode allows carers to carry out personal care in safety and comfort; minimising the risk of injury. Adjustment of the height and tilt is possible at the push of a button, ensuring an effortless working environment.

The **Ocean E-VIP** can be accessorised to suit a wide range of physical needs, making it a cost effective alternative to prescription based chairs. The **Ocean E-VIP** has a 150kg weight limit.



### Ocean E-VIP

- Seat height can be adjusted between 500 - 1050mm to allow easy access and effortless transfers.
- Seat tilt can be tilted forwards and backwards (tilt angle  $-5^{\circ}$  to  $35^{\circ}$ ) to achieve a comfortable seating position. The  $-5^{\circ}$  position offers easier transfer into the chair.
- Particularly stable even at the maximum tilt angle and maximum seat height.
- Clearly laid out hand control for simple adjustment of the seat height and tilt (4 buttons only).
- Fully sealed components: Safe for use in the shower.
- Can be cleaned in automatic washers at temperatures up to  $85^{\circ}$  C.
- 3 year warranty.



Comfortable powered recline



Optimal eye contact

### Standard features



#### Directional lock

Rear castor is equipped with a directional lock for straight travel to reduce back strain.



#### Headrest

Comfortable headrest which can be adjusted horizontally and vertically.



#### Tension adjustable back

Machine washable adjustable backrest, can be adjusted to the client.



#### Toilet pan

With handle and cover.

Aquatec **Ocean E-VIP options and accessories**



**Standard Soft Seat**

With key shaped cut out for hygiene access.  
**Supplied as standard.**



**Soft seat insert**

For soft seat with key-shaped cut out when commode function is not required.  
**Supplied as standard.**



**Small soft seat**

With small opening suitable for smaller and young adults.



**Splash guard / pommel**

Simply snap-on. (Fits soft seat and key-shaped cut out).



**Universal soft seat**

Closed seat surface. Smooth and comfortable.



**Safety bar**

For additional safety. Client can grip bar for additional re-assurance.



**Chest-belt & Pelvic Belt**

To ensure necessary safety during the positioning process. Chest belt standard. Pelvic belt optional.



**Armrest locking device**

Locks armrest into upright or lowered position.



**Lateral support**

This support reduces the backrest width and can help position the trunk.



**Soft backrest cushion**

For additional stability positions. Will benefit small individuals and young adults.



**Wider XL armrests**

Increases width between armrests by 40mm each side. Does not increase the weight limit.



**Calf & amputee support**

Padded. Both height and angle adjustable.



**Calf support with footrest**

Increased support helps distribute load to avoid pressure on lower limbs.



**Shortened leg rests**

Can be used to support abnormal posture or smaller/young adults.  
**Can be used across all models.**



**Ergonomic headrest holder**

Allows various adjustment positions to manage the clients head position.



**Head/Neckrest**

For ergonomic headrest holder. Offers exceptional management of head position. (Neckrest shown)

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Seat height	Seat width	Total width	Total depth	Total height	Width between armrests	Tilt angle seat	Max. user weight
<b>Ocean E-Vip</b>	500 - 1050 mm	450 mm	730 mm	1010 mm (without footrests) 1083 mm (with footrests)	1070 - 1620 mm (without headrest)	470 mm	0° - 35°	150 kg



**Ocean E-Vip** 39 kg

## Warranty

3 years warranty

\* A full risk assessment should be conducted prior to use with one carer.

## Did you know?

Invacare's Business Development Managers can provide support with **Ocean** Shower commode assessments. This service is FREE of charge. Contact your Invacare Business Development Manager to find out more.

E-mail [uk@invacare.com](mailto:uk@invacare.com)

## Have you seen the *Aquatec Ocean* support pack?

This resource is designed to provide additional straightforward, easy to understand information on the entire **Ocean** shower chair range. The pack includes a number of resources such as:

- Quick guide (A5 size)
- Ocean poster
- Ocean brochure
- Order form

To receive your **FREE** *Aquatec Ocean* support pack contact Invacare today by emailing [uk@invacare.com](mailto:uk@invacare.com)



Safe  
Patient  
Handling

The **Ocean E-VIP** is part of the Safe Patient Handling range.

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

