



The sturdy bath board

The *Aquatec Bali* serves for a safe seat on all standard bathtubs. It allows the user to safely and quickly wash with the help of the adjustable shower head whenever it is preferred to a bath.





Features and Options

Sturdy, stable and secure

Load capacity of 175 kg!

Safe and comfortable use

Rounded edges facilitate secure transfer. The non-slip structured surface provides safe sitting.



Adaptable

Secure fixation on different sizes of bathtubs with adjustable stoppers.



Stable handgrip

Offers additional safety.



Compartment for soap

Convenient and practical.

Technical data

Aquatec Bali	69,5 cm	30 cm	47 - 61 cm	22 cm	1,8 kg	175 kg

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
h@europa@invacare.com
www.invacare.eu.com

© 2015 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Bali - EU - 08/2015

Frame colour



White-Blue



Yes, you can.®