



Trusted Solutions,  
Passionate People

Power Wheelchair  
IRELAND Prescription Form  
PRONTOM41PP010117

Tel: (00353) 1 8107084 Fax: (00353) 1 8107085  
Email: ordersireland@invacare.com  
Online Spares available at www.invacare.ie

JANUARY 2017

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



# Pronto M41 Ireland

Max user weight 136kg (21.4 stone)



## Powered Positioning

### Standard Features



- Max user weight 136kg (21.4 stone)
- Compact, manoeuvrable, **MWD Pronto M41 Base**
- Tension Adjustable Backrest
- Telescopic adjustable Seat Widths from 380mm - 580mm (15" - 23") additional 50mm (2") by arms
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Detachable T-Armrests, height, width and depth adjustable with standard armpad
- Powered seat tilt (0° - 30°) and manually settable backrest recline (0° - 30°)
- Swing away and detachable legrest hangers with fixed footplate
- ACS DX2 Expandable Control System with REM A Joystick with standard bracket
- 255mm x 80mm (10" x 3") Solid Insert/Puncture Proof Grey Drive Wheels
- Candy Red, Deep Blue Shroud Colour
- 6km/h (4mph) In-line Motors and 32Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Toolkit and 8A (200 - 250V) Charger



**Standard Base Price : €4,143**

\* Configurations may vary from the image shown

X		Please make your selection with a cross in the box		◆ Standard feature	○ Free of Charge
Code		Description			Ireland
AQD0059		Pronto M41 "Modulite"		◆	
AQX0059		Pronto M41 "Modulite" Special Basic Model-PLEASE COMPLETE DETAILS AT THE END OF THE FORM			POA
<b>SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM</b>					
AQD2100	X	Modulite Seating Unit			€ 4,143
<b>SEAT FRAME AND SEAT PLATE</b>					
AQD2111	X	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)		◆	

**MODULITE SEAT SIZE SETTING CHART:**

			SEAT WIDTHS				
			330mm* (13")	AQD2121 380mm (15")	AQD2122 430mm (17")	AQD2123 480mm (19")	AQD2124 530mm (21")
Seat Depth AQD2129	Standard Seat Rail 410 - 510mm	AQD2132 (16") 410mm	○	○	○	○	○
		AQD2134 (18") 460mm	○	○	○	○	○
		AQD2136 (20") 510mm	○	○	○	○	○
Seat Depth AQD2130	Extended Seat Rail 510 - 560mm	AQD2136 (20") 510mm	○	○	○	○	○
		AQD2138 (22") 560mm	○	○	○	○	○

*Invacare* **Modulite**<sup>™</sup>  
An evolution in seating designed around you.



\*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

1. Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)

2. No seat cushion available with 560mm seat depth (2138)

3. Seat width 530mm (2124) not possible in combination with following armrests (2192/4)

4. Seat width 560mm (2138) When in combination with Flex Tension Adjustable Backrest the seat depth is reduced by 10mm

**SEAT CUSHION (sizes according to seat frame settings)**

AQD2142	Matrix Flo-tech Lite - pressure relieving seat cushion	€ 114
AQD2143	Matrix Flo-tech Contour Visco 430x430 mm	€ 114
AQD2147	Matrix Libra Cushion, 90mm thick	Not available with Seat Width 530mm with 410mm seat Depth, seat depth 560mm only available with 530mm seat width € 426
AQD2140	Without seat cushion	○

**BACKREST SUPPORT**

AQD2162	Backposts 540mm (tall) (must be selected for Flex3 adjustable backrest plates)	◆
AQD2163	Backposts 480mm (short)	Not available with Flex3 Adjustable backrest ○

**PUSH HANDLES FOR BACKPOSTS**

AQD2107	Backposts WITH push handles	◆
AQD2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. If you require a specific length of push handle please note in the specials box at the end of the form.	Not available with attendant controls, dual controls, Matrix backrest or tension adjustable backrest in conjunction with a headrest. € 44

**BACKREST UPHOLSTERY / CUSHION**

AQD2171	Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear	Not Available with Laterals, black fabric upholstery only, Headrest option AQD1496 ◆
AQD2179	Flex3 Tension Adjustable (long = 540mm / short = 480mm), available in Black Dartex only	Not available with Flip up or Following armrests, Chin Remote and Soft Armpads € 165
AQD2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)	€ 245
AQD2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	€ 397
AQD2174	Matrix Posture Back, height 400mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, Headrest option AQD1496 € 543
AQD2175	Matrix Posture Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, Headrest option AQD1496 € 543
AQD2176	Matrix Elite Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, Headrest option AQD1496, seat width of 380mm (15") AQD2121 in combination with flip up or following armrests € 543
AQD2177	Matrix Elite TR Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, Headrest option AQD1496 € 543
AQD2178	Matrix Elite Deep Back, height 500mm (Mesh fabric)	Not available with Laterals or Headrest option AQD1496 € 602
AQD2170	Without Backrest Upholstery (Straight 1" diameter Back Posts/Back Frame only suitable for mounting specialist backrests)	Not available with laterals or Headrest option AQD1496 ○

**UPHOLSTERY COLOUR**

AQD2151	Black fabric (nylon 2-way-stretch material)	◆
AQD2152	Black Dartex	○
AQD2153	LK Version (black with orange seam)	○
AQD2154	Black Plush	○

<b>ARMREST (NOTE: HIP GUIDES ONLY POSSIBLE WITH FLIP UP AND FOLLOWING ARMREST)</b>			
AQD2191	Single post 'T-Armrest' removable and height adjustable (with stepless positioning relative to the backposts) Range of adjustment from seat plate 245mm to 315mm		◆
AQD2181	Extended height Single post 'T-Armrest': Range of adjustment from seat plate 295mm to 365mm		€ 44
AQD2183	In-rigged Armrests (reduces seat width by 25 mm each side)	Not Available with hip support	€ 26
AQD2193	Flip-up Armrest (height and depth adjustable)	Not Available with laterals or in combination with 530mm seat plate and Flex3 adjustable backrest plates.	€ 63
AQD2192	Following Armrest (height adjustable from 230 - 300 mm)	Not Available with laterals or in combination with 530mm seat plate and Flex3 adjustable backrest plates.	€ 181
AQD2194	Following Armrest (height adjustable from 300 - 360 mm)	Not Available with laterals or in combination with 530mm seat plate and Flex3 adjustable backrest plates.	€ 181
<b>ARMPADS</b>			
AQD2195	Short Armpad (325 mm)		◆
AQD2196	Long Armpad (375 mm)	Not Available with Seat depth 410mm (16")	○
AQD2197	Long Ergonomic Armpad (397 mm)	Not Available with Seat depth 410mm (16")	○
AQD2199	Soft Armpad (REA Azalea Wide Pad)	Not Available with Seat depth 410mm (16")	€ 93
<b>HIP SUPPORTS LH (for flip up and following armrests only)</b>			
AQD2184	Multi adjustable hip support LH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
AQD2186	Multi adjustable hip support LH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
AQD2188	Multi adjustable hip support LH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
<b>HIP SUPPORTS RH (for flip up and following armrests only)</b>			
AQD2185	Multi adjustable hip support RH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
AQD2187	Multi adjustable hip support RH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
AQD2189	Multi adjustable hip support RH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests	€ 247
<b>SEAT ACCESSORIES</b>			
<b>HEADREST/NECKREST (Dartex) - Adequate headrest option must be selected</b>			
AQD1350	Neckrest: height and depth adjustable (Azalea type)		€ 77
AQD1432	Headrest: height, depth and sideways adjustable (Azalea type)		€ 215
AQD1427	Headrest with cheek support: height and depth adjustable (Azalea type)		€ 215
AQD1496	All bracketry as option 1432 without headrest cushion		€ 80
AQD1495	Non Invacare Headrest - no mounting hardware (with Backbrace/pushbar dia 22mm in conjunction with telescopic seat frame/backposts)		○
<b>LATERAL SUPPORT</b>			
AQD1446	Lateral Support with extended arm, rigid pad (Dartex) (Azalea type) (Pair)		€ 160
AQD1445	Lateral Support with extended arm, multi-adjustable pad (Dartex)(Azalea type) (Pair)		€ 160

## SEAT LIFTER / SEAT TILT

		SEAT HEIGHT		
		Tilt Range	420mm (16")	Selection
AQD1093	Electric Tilt with COG Shift	0° - 30°	◆	
AQD1094	Seat Lifter with electric tilt (COG Shift), <b>Not possible with Rem A Control</b>	0° - 30°	€1,950	



## BACKREST RECLINER

		Recline Range		Selection
AQD0975	Manually settable recline	0° - 30°	◆	
AQD0977	Powered Recliner with quick release for easy folding of backrest	0° - 30°	€303	



## PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS


### CENTRE MOUNTED LEGRESTS (INCLUDING CALF PAD AND ANGLE ADJ FOOTPLATE 240mm depth)

#### FOOTREST HANGER TYPE (INCLUDING CALF PAD AND ANGLE ADJ FOOTPLATE)

AQD0556	Manual : Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 390 - 510mm	€ 291
AQD0589	LNX Powered : Centre mounted legrest, length compensating	€ 1,739

#### SIDE MOUNTED LEGRESTS

AQD0526	Standard hanger 80° Pair	Not Available with Angle adjustable footplate	◆
AQD0503	<b>VARI F</b> : Angle adjustable legrest hanger (swing away and detachable) (adjustable via tools 0° - 70°) : knee to heel length 290 - 460mm		○
AQD0533	<b>VARI A</b> : Manually angle adjustable legrest hanger (swing away and detachable) (adjustable via hand wheel / knob 0° - 70°) - complete with calf pad support : knee to heel length 290 - 460mm	Not Available with Fixed footplate	€ 121
AQD0532	<b>ADM Manual</b> : Manually angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (adjustable via quick release lever) - complete with calf pad support : knee to heel length 290 - 460mm	Not Available with Fixed footplate	€ 210
AQD0531	<b>ADE Powered</b> : Electrically angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (powered via joystick) - complete with calf pad support : knee to heel length 290 - 460mm	Not Available with Fixed footplate	€ 622









#### FOOTPLATE TYPE; PAIRS

AQD0512	Fixed Footplate	◆
AQD0504	Angle Adjustable Footplate	○
AQD0513	Angle and depth adjustable Footplate (Adapt)	€ 110

## REMOTE DX2 SYSTEM

### USER CONTROL + BRACKET

AQD1638	REM Remote Type A (enables operation up to 4 electric. adjustments)	◆
AQD1601	ACS REM 24SD	€ 259
AQD1614	G90A	€ 191
	Head-Remote, ultrasonic proximity sensor(only available as special) Please contact UK Specials Dept	POA
AQD1207	Standard Bracket RH	○
AQD1208	Standard Bracket LH	○
AQD1707	Swing Away + height adjust.RH High End	€ 137
AQD1708	Swing Away + height adjust. LH High End	€ 137

<b>ATTENDANT REMOTE + BRACKET</b>			
AQD1643		REM Type A (Up to 4 electric actuators)	€ 17
AQD1602		ACS REM 24SD	€ 205
AQD1616		G90A	€ 205
AQD1207		Standard Bracket RH	○
AQD1208		Standard Bracket LH	○
<b>CHIN CONTROL + BRACKET</b>			
AQD1719		Chin Remote with standard Joystick (Always to be mounted on the Master Remote's side) (for Modulite only)	(Only compatible with Option1601 ) € 2,184
AQB1289		Swing Away Bracket RH, side manual	Only available with RH user bracket ○
AQD1288		Swing Away Bracket LH, side manual	Only available with LH user bracket ○
AQD1287		Swing Away Bracket RH, powered	Only available with RH user bracket € 1,472
AQD1286		Swing Away Bracket LH, powered	Only available with LH user bracket € 1,472
<b>DUAL CONTROL + BRACKET (can only be selected as an addition to user remote)</b>			
AQD1720		Compact Dual Control	€ 471
AQD1207		Standard Bracket RH	Only available with LH user bracket ○
AQD1208		Standard Bracket LH	Only available with RH user bracket ○
<b>INVACARE IDC DUAL CONTROL + BRACKET (Only available with ACS)</b>			
AQD1645		Invacare IDC Dual control	€ 1,256
AQD1221		Detachable Bracket RH	○
AQD1222		Detachable Bracket LH	○
<div style="display: flex; justify-content: space-around; text-align: center;"> <div> <p><b>REM Type A</b></p>  </div> <div> <p><b>REM 24SD</b></p>  </div> <div> <p><b>G90 A</b></p>  </div> <div> <p><b>Chin control</b></p>  </div> <div> <p><b>Dual control</b></p>  </div> <div> <p><b>IDC Control</b></p>  </div> </div>			
<b>DRIVE UNIT OPTIONS</b>			
<b>FRAME COLOUR (COLOURED PARTS IS A DETACHABLE AND CHANGEABLE PLASTIC SHROUD)</b>			
AQD0227		Candy Red	◆
AQD0283		Deep Blue	○
<b>CENTRE WHEELS (DRIVING WHEELS)</b>			
AQD0781	X	Centre Wheel 255mm x 80mm (10"-3") Solid Insert/Puncture Proof Grey	◆
<b>PERFORMANCE OPTIONS</b>			
<b>MOTOR TYPE</b>			
AQD1282	X	6 km/h (4mph)	◆
<b>BATTERY</b>			
AQD1261	X	Battery 32Ah AGM pair	◆
<b>CHARGER</b>			
AQD1343	X	Charger 8A (200 - 250V)	◆
<b>SAFETY ACCESSORIES</b>			
<b>TRANSPORTATION KIT</b>			
AQD1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) <b>(Adequate headrest is compulsory)</b>	◆
<b>FIXATION BELT</b>			
AQD1329		Lapbelt (adjustable both sides) fitted to base of backrest	○
AQD2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	Not Available with seat depth 410mm (16") and T-Armrest ◆
<b>ACCESSORIES</b>			
<b>CANE HOLDER</b>			
AQD1312		RH	€ 37
AQD1363		LH	€ 37
<b>TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)</b>			
AQD1404		Tray Swing Away RH	€ 333
AQD1405		Tray Swing Away LH	€ 333
<b>TOOLKIT</b>			
AQD1333	X	Toolkit (incl bit holder, bits and fork spanners)	◆
<b>TOTAL PRICE</b>			





Need a different option not on the form?

Take a look at our Powerchair Alternative Options form, available as standard.

Powerchair Alternative Option required:



**Special Instructions - Please supply the following information:**

**Name:**

**Email:**

**Phone:**

**Fax:**

**Description of special request (please include relevant option code or part number):**