



Invacare® | ***Matrix® Flo-tech™ Lite Visco***

Low profile visco elastic memory foam cushion offering exceptional comfort



Trusted Solutions,
Passionate People



The *Invacare* **Matrx Flo-tech Lite Visco** cushion is made of high-density visco elastic foam, allowing it to conform to the user's shape. This affords exceptional levels of comfort and increases and redistributes pressure. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.



Features and Options



Low memory foam designed to maximise body contact area

The visco-elastic foam conforms to the individual's body, maximising the contact area, optimising pressure distribution and delivering exceptional levels of comfort.

Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle for easy transportation.

Available in 61 sizes

The *Invacare* **Matrx Flo-tech Lite Visco** is available in a range of 61 sizes, from 250 x 250 to 480 x 530 mm. (See prescription chart for details).

Three foam densities available







Available in three densities to accommodate a range of client weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 100 kg).

Contoured seating surface with pommel and side support

Offers mild lateral support and pelvic stability to enhance the individual's posture.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see the *Invacare* website.

| Various sizes | Height | Total product weight | Max. user weight ¹ | Max. cleaning temperature | Warranty | |
|---|---|---|---|--|---|---|
|  |  |  |  |  |  | Manufactured to European Standards: Interior Foam: BS EN 1021-1, BS EN 1021-2, BS 5852: 1982 Crib 5 |
| Matrx Flo-tech Lite Visco A range of 61 sizes | 78 mm | 1.26 kg | 100 kg (16 stone) | 80° | 3 years | |

Disclaimer

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

Colours



Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Limited
Pencoed Technology Park
Pencoed - Bridgend CF35
5AQ - United Kingdom

Tel. +353 1 8107084
Fax +353 1 8107085
Email: ireland@invacare.com
www.invacare.ie

Invacare Ireland Ltd
Unit 5 - Seatown
Business Campus
Seatown Road
SWORDS - County
Dublin - Ireland

