



Invacare® | **Matrix® Flo-tech™ Contour Lo-Back**

Allows the individual to sit back further into their seat

**matrix**  
SEATING SERIES

Trusted Solutions,  
Passionate People

**INVACARE**  
Yes, you can.®

The *Invacare **Matrix Flo-tech Contour Lo-Back Cushion*** has been designed with a lower rear contour allowing the individual to sit further back into seat. The additional contouring and support under the thighs stabilizes the individual further, thus reducing the damaging effects of shear.

### Comparison



*Matrix Flo-tech Contour*



*Matrix Flo-tech Contour Lo-Back*

### Contour Lo-Back Cushion optimizes support surface contact area

Comparison with standard wheelchair cushion, showing improved posture through use of the *Matrix Flo-tech Contour Lo-Back Cushion*.

## Features and Options



### Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle for transportation.



### Available in 83 sizes

The *Invacare **Matrix Flo-tech Contour Lo-Back Cushion*** is available in a range of 83 sizes, from 360 x 330 to 600 x 550 mm. See size chart for details.



### Contoured seating

Offers improved stability and postural support. Castellations reduce shear and friction forces on tissue.



### Ultra-fresh™ added to all cushions for antimicrobial and odour protection

Antimicrobial treatment has been added to all ***Flo-tech*** foam to help make it resistant to the growth of odour causing bacteria and helps provide longer lasting freshness.

### Moulded Foam Cushion

Moulding allows ***Matrix Flo-Tech*** cushions to provide a three-dimensional form, which cannot be emulated in cut foam. Moulded foam is highly durable, increasing the product's longevity.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see the [Invacare website](#).

Various sizes	Height	Total product weight	Max. user weight <sup>1</sup>	Max. cleaning temperature	Warranty
A range of 83 sizes	103 mm	1.38 kg	152 kg (24 stone)	80°	3 years


Manufactured to European Standards: BS EN 1021-1 Cig, BS EN 1021-2 Match and BS 7176 (Crib 5)

#### Disclaimer


1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

## Colours

Black Smooth



Black wipedown



Black spacer fabric



Tel: +44 1 656 776 222  
 Fax: +44 1 656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Limited  
 Pencoed Technology Park  
 Pencoed - Bridgend CF35  
 5AQ - United Kingdom

Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [ireland@invacare.com](mailto:ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

Invacare Ireland Ltd  
 Unit 5 - Seatown  
 Business Campus  
 Seatown Road  
 SWORDS - County  
 Dublin - Ireland



**Yes, you can.®**

