

# ***Flo-tech<sup>®</sup> Image***



Slimline cushion with fluid sac



Yes, you can.

The *Invacare Flo-tech Image* cushion uses a combination of a slim, shaped foam cushion base and a fluid sac covering the sacral area, increasing pressure distribution and immersion of the ischial tuberosities. This design delivers a high degree of protection to users who are at High Risk of developing pressure ulcers. This slim cushion format is ideal for both chair and wheelchair use.

## Features and Options



### Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.

### Available in 23 sizes

The *Invacare Flo-tech Image* is available in a range of 23 sizes, from 360 x 330 to 460 x 410 mm. (See size chart for details).

### Integral Sag Compensator option

Option of sag compensator to accommodate the sag in wheelchair seats, and ensure the cushion is correctly presented.

### Ultra-fresh™ added to all cushions for antimicrobial and odour protection

Antimicrobial treatment that has been added to all **Flo-tech** foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.

### Slimline contoured foam and gel combination cushion

Some users prefer a slimline cushion, finding it more compatible with their daily living activities. The *Invacare Flo-tech Image* combines the benefit of a contoured foam upper surface with a large fluid sac covering the ischial and sacral area, offering increased support and reducing the peak pressures.

### Contoured seating surface with gel sac

Offers improved stability and postural support; moulded foam castellations and fluid combine effectively to reduce shear and friction forces on tissue.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Various sizes	Height	Total product weight	Max. user weight <sup>1</sup>	Max. cleaning temperature
<b>Flo-tech Image</b> A range of 50 sizes	78 mm	1.82 kg	140 kg	80°

Manufactured to European Standards:  
BS EN 1021-1 Cig, BS EN 1021-2 Match  
and BS 7176 (Crib 5)

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

## Colours

Black Smooth	Black wipedown	Black spacer fabric	Baltic Blue

T: +44 1656 776222  
uk@invacare.com  
www.invacare.co.uk  
www.invacare.ie

Invacare UK & Ireland  
Unit 4 Pencoed  
Technology Park  
Brigend CF35 5AQ



**Yes, you can.™**