

Essential Visco

Quality Foam Mattress



A cost effective solution for comfort
and pressure redistribution



The **Essential Visco** mattress is a pressure redistributing support surface for the home care, residential and long-term care environment. This mattress features a high quality multi-stretch vapour-permeable, smooth and water resistant, polyurethane, fabric cover that helps to reduce shear and friction forces on the delicate skin of the patient, in addition to helping meet modern hygiene standards.

- ▶ Visco elastic foam top layer
- ▶ Sturdy foam base for additional support
- ▶ Multi-stretch polyurethane fabric cover
- ▶ Zipped around 3 sides for ease of inspection and laundering



Visco Elastic Foam Top layer

- ▶ Conforms to the body shape offering maximum support and pressure relief.

Cover

- ▶ The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability, which helps to reduce shear and friction forces. The cover also has a white inner liner allowing easy inspection for strikethrough or other cover damage.

Design

- ▶ The **Essential Visco** mattress has an innovative and modern cover design with the simple wave giving the mattress a fresh and current feel.

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Size range ▼



Max. user weight ▼



Max. cleaning temperature ▼

Nominal density range (kg/m³) ▼

Article no. ▼

	Size range ▼	Max. user weight ▼	Max. cleaning temperature ▼	Nominal density range (kg/m ³) ▼	Article no. ▼
Essential Visco (150 mm height)	2000 x 880 x 150 mm	140 kg	95°	32 - 34, Visco layer 48 - 52	1588533
Essential Visco (120 mm height)	2000 x 880 x 120 mm	110 kg	95°	32 - 34 Visco layer 48 - 52	1588538

Mattresses is manufactured to:
EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

1. *Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).*

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie



Yes, you can.®

