



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 2NG Stock Spec 1

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Buckle belt (adjustable both sides)
- Anti-tippers (24" wheels only)
- Anthracite Black
- Choice of 5 seat widths 15" - 19"
- Depth adjustable 16" to 17"
- Black padded nylon seat upholstery
- Half folding 7° backrest 17" (430mm)
- Swing in/out 80° detachable hanger with adjustable heel strap
- Removable, lightweight and height adjustable full arm (200mm > 260mm)
- Black 24" x 1.3/8" quick release rear wheel solid tyre or
- 12" X 2" fixed axle, black solid tyre wheel forward position
- Black 8" x 1.1/4" (200mm x 25mm) solid front castor aluminium fork
- Push to lock brakes



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
A2NG SELF PROPEL		
BQ1529465-AF	A2NG 15" W X 16" to 17" D	€331
BQ1529466-AF	A2NG 16" W X 16" to 17" D	€331
BQ1529469-AF	A2NG 17" W X 16" to 17" D	€331
BQ1529470-AF	A2NG 18" W X 16" to 17" D	€331
BQ1529485-AF	A2NG 19" W X 16" to 17" D	€331

A2NG TRANSIT

BQ1529486- AF		A2NG 15" W X 16" to 17" D	€331
BQ1529487- AF		A2NG 16" W X 16" to 17" D	€331
BQ1529490- AF		A2NG 17" W X 16" to 17" D	€331
BQ1529491- AF		A2NG 18" W X 16" to 17" D	€331
BQ1529492- AF		A2NG 19" W X 16" to 17" D	€331

TOTAL PRICE	
--------------------	--