



Yes, you can:

Manual Invacare Ireland Prescription Form

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3 Junior Basic

Standard Features

- Max user weight: 80kg (12.5 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 3 seat widths 12" - 14" (305mm to 355mm)
- Choice of 2 seat depths (or depth adjustable upholstery)
- Choice of 3 seat heights 16" - 18" (410mm to 460mm)
- Black nylon upholstery
- Folding backrest set at 16" (400mm)
- Width and height adjustable backrest upholstery
- 80° detachable legrest hangers
- Flip up, angle and depth adjustable footplates
- Removable and height adjustable desk armrests
- Adjustable heel straps
- Choice of 7 colours (with 8 special colour options)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt
- User tool kit

Standard Base price: €865



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CNB0000	ACTION ³ JUNIOR	○
CNBSPEC	SPECIAL	POA
PROPULSION*		
CNB0100	SELF PROPULSION (NC 1010, 1170)	€865
CNB0110	TRANSIT PROPULSION, WITH RH TIPPING AID (NC 1020, 1030, 1180, 1190, 1195, 1275, 1390, 1430, 1440)	€865

SEAT WIDTH*				
CNB0180		12" SW	(NC 0250, 0252, 0520, 0570, 1190, 1195)	○
CNB0190		13" SW	(NC 0240, 0251, 0530, 0580)	○
CNB0200		14" SW	(NC 0240, 0251, 0530, 0580)	◆
SEAT DEPTH*				
CNB0240		12" SD		○
CNB0250		14" SD		○
CNB0251		ADJ SD 10" - 12"	(NC 1364)	○
CNB0252		ADJ SD 12" - 14"		◆
SEAT HEIGHT*				
CNB0260		18" SH		○
CNB0270		17" SH		○
SEAT UPHOLSTERY*				
CNB0330		NYLON WITH HOOK AND LOOP STRIPS		◆
SEAT CUSHION*				
CNB33333		WITHOUT SEAT CUSHION		○
CNB0362		BASIC SEAT CUSHION (ESSENTIAL JUNIOR)		€132
BACKREST*				
CNB0370		FIXED 0° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM		○
CNB0380		FIXED 7° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM		○
CNB0390		RECLINING (RATCHET) 0/30° 5 HEIGHTS 13" - 17"	(NC 1190, 1195)	€109
CNB0395		RECLINING (GAS) 0/30° 16"	(NC 1190, 1195)	€180
CNB0400		FOLDING BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400 OR 17" / 430MM		◆
CNB0401		FOLIDNG BACKREST 20"		○
BACKREST UPHOLSTERY*				
CNB0465		WITHOUT BACK UPHOLSTERY	(NC 0855)	○
CNB0475		WIDTH AND HEIGHT ADJUSTABLE NYLON SLING TYPE		◆
CNB0480		CONTOUR (TENSION) ADJUSTABLE NYLON TYPE	(NC 0855)	€130
LEGREST HANGER RH*				
CNB0505		90° STRAIGHT RH REMOVABLE HANGER UL STD 125 TO 260 MM	(NC 0600)	€34
CNB0510		80° SHORT RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC 0562)	◆
CNB0512		80° LONG RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	(NC 0560)	○
CNB0520		ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM	(NC 0600, 0611, 0880, 0886)	€106
CNB0530		ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM	(NC 0600, 0611, 0880, 0886)	€106

LEGREST HANGER LH*			
CNB0555		90° STRAIGHT LH REMOVABLE HANGER UL STD 125 TO 260 MM (NC 0620)	€34
CNB0560		80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	◆
CNB0562		80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	○
CNB0570		ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM (NC 0611, 0620, 0880, 0886)	€106
CNB0580		ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM (NC 0611, 0620, 0880, 0886)	€106
FOOTREST RH*			
CNB0600		FLIP UP RH FOOTPLATE (NC 0875, 0886)	○
CNB0605		FLIP UP DEPTH AND ANGLE ADJ RH	◆
CNB0611		ONE PIECE FOOTBOARD (NC 0620, 0625, 0875)	○
FOOTREST LH*			
CNB0620		FLIP UP LH FOOTPLATE (NC 0875, 0886)	○
CNB0625		FLIP UP DEPTH AND ANGLE ADJ LH	◆
CNB66666		NO FOOTPLATE REQUIRED (PLEASE SELECT IN COMBINATION WITH OPTION 0611) (NC 0600, 0605, 0620, 0625)	○
ARMREST RH*			
CNB0640		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM (NC 0745, 0746, 0750, 0755)	◆
CNB0645		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM (NC 0740, 0746, 0750, 0755)	○
CNB0646		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM (NC 0745, 0750, 0755)	○
CNB0650		MUDGUARD RH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0740, 0745, 0746, 0755, 0910, 1201, 1471, 1472)	€35
CNB0655		SKIRTGUARD RH JUNIOR (NC 0740, 0745, 0746, 0750, 0910, 1471, 1472)	○
ARMREST LH*			
CNB0740		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM	◆
CNB0745		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM	○
CNB0746		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM	○
CNB0750		MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0910, 1201, 1471, 1472)	€35
CNB0755		SKIRTGUARD LH JUNIOR (NC 0910, 1471, 1472)	○

OTHER SEATING OPTIONS				
CNB0855		BACKREST EXTENSION	(NC 0860, 0865)	€85
CNB0860		HEADREST STANDARD MULTI - ADJUSTABLE	(NC 0865, 0902)	€130
CNB0865		HEADREST ROUND MULTI - ADJUSTABLE	(NC 0902)	€182
CNB0880		ADJUSTABLE CALF STRAP		€13
CNB0881		FRONT STABILISER KIT + 100 MM		○
CNB0875		BODYPOINT ANKLE HUGGERS (PAIR) SIZE S - UNMOUNTED		€189
CNB0886		FOOTPLATE ADJUSTMENT OPTIMIZER KIT		€24
CNB0902		WIDTH ADJUSTABLE BACK BRACE		€30
CNB0910		WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1		○
COLOUR*				
CNB0952		PEARL WHITE		○
CNB0953		MATT BLACK		○
CNB0954		MATT RED		○
CNB0955		INTENSE YELLOW		○
CNB0956		POLAR BLUE		○
CNB0957		CANDY PINK		○
CNC0959		BLUE		○
CNB0991		SPECIAL COLOUR (ON DEMAND ONLY)		€164
SWO0001		SLATE GREY - 3 WEEK LEAD TIME		€164
SWO0002		ELECTRIC GREEN - 3 WEEK LEAD TIME		€164
SWO0003		HYPNOTIC PURPLE - 3 WEEK LEAD TIME		€164
SWO0004		ANTHRACITE BLACK - 3 WEEK LEAD TIME		€164
SWO0005		DARK RED - 3 WEEK LEAD TIME		€164
SWO0006		PEARL GREY - 3 WEEK LEAD TIME		€164
SWO0007		VIBRANT ORANGE - 3 WEEK LEAD TIME		€164
FRONT CASTORS*				
CNB0999		5" X 1" BLACK SOLID TYRE		◆
CNB1001		5" X 1" FLASHING LIGHTS SOLID TYRE		€37
CNB1007		5" X 1.3/4" BLACK SOLID TYRE		○
REAR WHEEL SIZE*				
CNB1010		12"	(NC 1045, 1130, 1140, 1141, 1180, 1253, 1254, 1255, 1256, 1380, 1390)	○
CNB1020		20"	(NC 1055, 1050, 1170)	◆
CNB1030		22"	(NC 1045, 1055, 1170)	○
REAR WHEELS*				
CNB1045		REINFORCED STANDARD SELF-PROPEL REAR WHEEL 6 POINT FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)	(NC 1280)	◆
CNB1050		DRUM BRAKE 4 POINT FIXING HANDRIM (22/24" COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES / 12" BLACK PLASTIC HUB AND RIM) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)	(NC 1270, 1275)	○
CNB1055		TRANSIT BLACK PLASTIC	(NC 1280)	○
TYRES*				
CNB1120		SOLID BLACK		◆
CNB1130		PNEUMATIC PROFILE BLACK		○

HANDRIM - A SELECTION MUST BE MADE FOR SELF-PROPULSION			
CNB1140	<input type="checkbox"/>	ALUMINIUM PRE-ANODIZED HANDRIM	○
CNB1141	<input type="checkbox"/>	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	€88
REAR WHEEL AXLE*			
CNB1170	<input type="checkbox"/>	FIXED AXLE	◆
CNB1180	<input type="checkbox"/>	Q/R AXLE	○
REAR WHEEL POSITION*			
CNB1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CNB1195	<input type="checkbox"/>	MEDIUM REAR WHEEL POSITION	○
CNB1200	<input type="checkbox"/>	PASSIVE REAR WHEEL POSITION	○
CNB1201	<input type="checkbox"/>	DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	€83
CAMBER ANGLE*			
CNB1205	<input type="checkbox"/>	For Transit and Drum brake	○
CNB1210	<input type="checkbox"/>	2° CAMBER (NC 1280)	○
CNB1211	<input type="checkbox"/>	3° CAMBER STANDARD (NC 1280)	◆
CNB1212	<input type="checkbox"/>	4° CAMBER (NC 1280)	○
CHASSIS OPTION			
CNB1253	<input type="checkbox"/>	SPOKE GUARDS CLEAR	€93
CNB1254	<input type="checkbox"/>	SPOKE GUARDS PIXEL	€113
CNB1255	<input type="checkbox"/>	SPOKE GUARDS SNOWFLAKES	€113
CNB1256	<input type="checkbox"/>	SPOKE GUARDS COLOR ME	€113
BRAKES*			
CNB1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CNB1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CNB1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES	€179
TRANSPORTATION KIT*			
CNB1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CNC1362	<input type="checkbox"/>	BODYPOINT POSTURE BELT 2 POINTS, SMALL 110 mm LENGTH – 25 mm WEBBING	○
CNC1364	<input type="checkbox"/>	BODYPOINT POSTURE BELT 2 POINTS, MEDIUM, 140 mm LENGTH – 40 mm WEBBING	◆
CNB1370	<input type="checkbox"/>	REAR ANTI-TIPPERS WITH TIPPING AID (NC 1430, 1440)	€41
CNB1380	<input type="checkbox"/>	REFLECTORS KIT	€10
CNB1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	€16
OTHER OPTIONS			
CNB1430	<input type="checkbox"/>	RH TIPPING AID	€16
CNB1440	<input type="checkbox"/>	LH TIPPING AID	€16
CNB1471	<input type="checkbox"/>	RH SWING AWAY ACRYLIC TRAY	€235
CNB1472	<input type="checkbox"/>	LH SWING AWAY ACRYLIC TRAY	€235

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CNB0240	12" SD
CNB0250	14" SD
CNB0251	ADJ SD 10" - 12"
CNB0252	ADJ SD 12" - 14"
CNB0520	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM
CNB0530	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM
CNB0560	80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY
CNB0562	80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM
CNB0570	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM
CNB0580	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM
CNB0600	FLIP UP RH FOOTPLATE
CNB0605	FLIP UP DEPTH AND ANGLE ADJ RH
CNB0611	ONE PIECE FOOTBOARD
CNB0620	FLIP UP LH FOOTPLATE
CNB0625	FLIP UP DEPTH AND ANGLE ADJ LH
CNB0740	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM
CNB0745	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM
CNB0746	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM
CNB0750	MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM)
CNB0755	SKIRTGUARD LH JUNIOR
CNB0855	BACKREST EXTENSION
CNB0860	HEADREST STANDARD MULTI - ADJUSTABLE
CNB0865	HEADREST ROUND MULTI - ADJUSTABLE
CNB0880	ADJUSTABLE CALF STRAP
CNB0875	BODYPOINT ANKLE HUGGERS (PAIR) SIZE S - UNMOUNTED
CNB0886	FOOTPLATE ADJUSTMENT OPTIMIZER KIT
CNB0902	WIDTH ADJUSTABLE BACK BRACE
CNB0910	WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1
CNB1010	12"
CNB1020	20"
CNB1030	22"
CNB1045	REINFORCED STANDARD SELF-PROPEL REAR WHEEL 6 POINT FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CNB1050	DRUM BRAKE 4 POINT FIXING HANDRIM (22/24" COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES / 12" BLACK PLASTIC HUB AND RIM) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CNB1055	TRANSIT BLACK PLASTIC
CNB1130	PNEUMATIC PROFILE BLACK
CNB1140	ALUMINIUM PRE-ANODIZED HANDRIM
CNB1141	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CNB1170	FIXED AXLE
CNB1180	Q/R AXLE
CNB1190	STANDARD REAR WHEELS POSITION
CNB1195	MEDIUM REAR WHEEL POSITION

CNB1201	DUAL AMPUTEE REAR WHEEL POSITION
CNB1253	SPOKE GUARDS CLEAR
CNB1254	SPOKE GUARDS PIXEL
CNB1255	SPOKE GUARDS SNOWFLAKES
CNB1256	SPOKE GUARDS COLOR ME
CNB1270	PUSH TO LOCK BRAKES
CNB1275	PULL TO LOCK BRAKES
CNB1280	ATTENDANT HUB BRAKES
CNC1364	BODYPOINT POSTURE BELT 2 POINTS, MEDIUM, 140 mm LENGTH – 40 mm WEBBING
CNB1380	REFLECTORS KIT
CNB1390	BRAKE EXTENSION PAIR
CNB1430	RH TIPPING AID
CNB1440	LH TIPPING AID
CNB1471	RH SWING AWAY ACRYLIC TRAY
CNB1472	LH SWING AWAY ACRYLIC TRAY



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):