



**Yes, you can:**

# Manual Invacare Ireland Prescription Form

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 4NG Transit

### Standard Features

- Attendant Propulsion
- Choice of 5 seat widths (15" - 19") 380mm to 480mm
- Choice of 3 seat depths (16" - 20") 400mm to 500mm
- Choice of 3 seat heights (18" - 20") 460mm to 510mm
- Black nylon seat upholstery
- Fixed 7° 380mm backrest
- Standard padded lumbar support nylon sling type backrest upholstery
- 80° Swing in/Swing out legrest hangers
- Flip up footplates
- Swing away detachable desk armrests
- Choice of 8 colours (with 7 special colour options)
- Push to lock brakes
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

**Standard Base price: €942**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
CKF0000	ACTION 4NG	○
CKFSPEC	SPECIAL	○
<b>PROPULSION*</b>		
CKF0110	TRANSIT PROPULSION	(NC 1111) €942

SEAT WIDTH*			
CKF0120		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm) (NC 0656)	○
CKF0130		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CKF0140		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CKF0150		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	◆
CKF0160		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CKF0180		16" (400mm) SD (NC 0702, 0703)	○
CKF0190		18" (450mm) SD (NC 0701, 0703)	◆
CKF0195		20" (500mm) SD (NC 0280, 0295, 0659, 0661, 0681, 0691, 0701, 0702)	○
CKF0200		16" - 18" ADJ SD (NC 0280, 0295, 0703)	€25
SEAT HEIGHT*			
CKF0210		20" (510mm) SH	○
CKF0220		19" (485mm) SH	◆
CKF0230		18" (460mm) SH	○
SEAT UPHOLSTERY*			
CKF0280		NYLON W/O HOOK AND LOOP STRIPS	○
CKF0290		NYLON WITH HOOK AND LOOP STRIPS	◆
CKF0295		PADDED NYLON	€25
BACKREST - Anti-tippers required with reclining backrests (0320, 0330 and 0334)*			
CKF0300		FIXED 7° BACKREST HT ADJ 15" - 18" (NC 0364, 0367, 0368, 0375, 0537, 0538, 0539, 0540, 0597, 0598, 0599, 0600, 0716, 0717, 0718, 0780)	◆
CKF0310		FIXED 7° BACKREST HT ADJ 15" - 18" HT ADJ P/HANDLES (NC 0364, 0367, 0368, 0375, 0537, 0538, 0539, 0540, 0597, 0598, 0599, 0600, 0716, 0717, 0718, 0780)	€105
CKF0311		FIXED 7° BACKREST 2 HEIGHTS 16" OR 17" (NC 0360, 0367, 0782)	○
CKF0312		FIXED 7° BACKREST 20" (NC 0360, 0367, 0782)	○
CKF0314		FIXED 7° BACKREST 17" ADJ P/HANDLES 6" (NC 0360, 0367, 0782)	€50
CKF0318		ANGLE ADJUSTABLE (-15° TO +15°) (NC 0350, 0360, 0364, 0368, 0375, 0716, 0717, 0718, 0721)	€101
CKF0320		RECLINING (RATCHET) 0/30° HEIGHTS 15"/16"/17" (NC 0360, 0367, 0780, 0782)	€152
CKF0334		RECLINING (GAS) 0/30° HEIGHTS 18"/19"/20" (NC 0367, 0780, 0782)	€226
CKF0341		FOLDING 7° BACKREST 17" (NC 0360, 0367, 0782)	○
CKF0342		FOLDING 7° BACKREST 20" (NC 0360, 0367, 0375, 0782)	○
BACKREST UPHOLSTERY*			
CKF0350		WITHOUT BACK UPHOLSTERY (NC 0721)	○
CKF0360		STANDARD PADDED NYLON SLING TYPE (JEMIMA) - AVAILABLE WITH BACKREST OPTIONS 0300 AND 0310 ONLY (NC 0721)	○
CKF0364		LUMBAR +1 STRAP	○
CKF0367		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL. BACK POCKET - AVAILABLE WITH BACKREST OPTION 0318 ONLY	€63
CKF0368		CONTOUR 3+1 STRAPS (NC 0780, 0782)	€37
CKF0375		CONTOUR 7+1 STRAPS (NC 0780, 0782)	€142
RH LEGREST HANGER*			
CKF0390		80° RH (NC 1111, 66666)	◆
CKF0400		ELEVATING LEGREST RH (NC 0491, 0760, 0770, 0771, 66666)	€150
CKF0420		AMPUTEES LEGREST RH (NC 0480, 0490, 0491, 0760, 0771, 1111)	€159

<b>LH LEGREST HANGER*</b>				
CKF0440		80° LH	(NC 2222)	◆
CKF0450		ELEVATING LEGREST LH	(NC 0491, 0760, 0770, 0771, 55555)	€150
CKF0470		AMPUTEE LEGREST LH	(NC 0491, 0500, 0510, 0760, 0771, 2222, 55555)	€159
<b>RH FOOTREST*</b>				
CKF0480		FLIP UP RH		◆
CKF0490		DEPTH AND ANGLE ADJUSTABLE RH		€50
SWO1111		FLIP UP DEPTH AND ANGLE ADJ RH (SEE FULL DESCRIPTION)		€50
CKF0491		ONE PIECE	(NC 0500, 0510, 0770, 0771)	€93
CKF66666		NO FOOTREST		○
<b>LH FOOTREST*</b>				
CKF0500		FLIP UP LH		◆
CKF0510		DEPTH AND ANGLE ADJUSTABLE LH		€50
SWO2222		FLIP UP DEPTH AND ANGLE ADJ LH (SEE FULL DESCRIPTION)		€50
CKF55555		NO FOOTREST		○
<b>RH ARMREST*</b>				
CKF0520		REMOVABLE RH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ.AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0593, 0596, 0597, 0598, 0599, 0600, 1270, 1271)	○
CKF0530		REMOVABLE RH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0593, 0596, 0597, 0598, 0599, 0600)	○
CKF0533		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0580, 0590, 0597, 0598, 0599, 0600, 1270, 1271)	€103
CKF0536		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0580, 0590, 0597, 0598, 0599, 0600, 1271)	€103
CKF0537		SWING AWAY DETACHABLE RH (N°2) DESK	(NC 0580, 0590, 0593, 0596, 0599, 0600, 1270, 1271)	◆
CKF0538		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH	(NC 0580, 0590, 0593, 0596, 0599, 0600)	○
CKF0539		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0580, 0590, 0593, 0596, 0597, 0598, 1270, 1271)	€67
CKF0540		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0580, 0590, 0593, 0596, 0597, 0598, 1271)	€71
CKF0555		REMOVABLE RH (N°4) ARMPAD HEMI MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED	(NC 1270, 1271)	€69

LH ARMREST*			
CKF0580		REMOVABLE LH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1270, 1272)	○
CKF0590		REMOVABLE LH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	○
CKF0593		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1270, 1272)	€103
CKF0596		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1272)	€103
CKF0597		SWING AWAY DETACHABLE LH (N°2) DESK (NC 1270, 1272)	◆
CKF0598		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH	○
CKF0599		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 1270, 1272)	€67
CKF0600		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 1272)	€71
CKF0615		REMOVABLE LH (N°4) HEMI MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED ARMPAD (NC 1270, 1272)	€69
SEAT CUSHION			
CKF0654		BASIC 50MM - NYLON AND ANT-SLIP COVER	€27
CKF0659		MATRX FLO-TECH LITE CUSHION - Low Risk (NC 0703)	€168
CKF0661		MATRX FLO-TECH IMAGE CUSHION - High Risk (NC 0703)	€382
CKF0681		MATRX FLO-TECH CONTOUR VISCO NG CUSHION - High Risk (NC 0703)	€236
CKF0691		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk (NC 0703)	€566
CKF0693		MATRX LIBRA CUSHION - Very High Risk	€645
SEAT CUSHION DEPTH			
CKF0701		16" SD	○
CKF0702		18" SD	○
CKF0703		20" SD (NC 0656)	○
POSITIONING ACCESSORIES FOR MATRX LIBRA & FLO-TECH			
CKF0706		TAPERED WEDGES (PAIR)	€28
CKF0707		ABDUCTOR	€58
CKF0708		LATERAL WEDGE	€46
CKF0709		ANTERIOR SEAT WEDGE	€43
OTHER SEATING OPTIONS			
CKF0656		SEAT RIGIDIZER	€79
CKF0716		NECKREST INCLUDING BACKBRACE (NC 0717, 0718, 0780, 0782)	€158
CKF0717		HEADREST INCLUDING BACKBRACE (NC 0718, 0780, 0782)	€190
CKF0718		HEADREST WITH CHEEK SUPPORT (NC 0780, 0782)	€217
CKF0721		LATERAL POSITIONING PADS PAIR	€72
SWO3333		MULT ADJUSTABLE LATERALS (NC 0300, 0310, 0314, 0318, 0320, 0334)	€340
CKF0760		CALF STRAP	○
CKF0770		LONGER FOOTREST TUBE (NC 0771)	€36
CKF0771		SHORTER FOOTREST TUBE	€48
CKF0780		BACKBRACE (NC 0782)	€24
CKF0782		FOLDABLE BACKBRACE	€45

COLOUR*			
CKF0801	<input type="checkbox"/>	SLATE GREY	○
CKF0804	<input type="checkbox"/>	PEARL WHITE	○
CKF0807	<input type="checkbox"/>	MATT BLACK	○
CKF0811	<input type="checkbox"/>	RED MATT	○
CKF0815	<input type="checkbox"/>	WATER BLUE	○
CKF0816	<input type="checkbox"/>	BLUE	○
CKF0817	<input type="checkbox"/>	DARK PURPLE	○
CKF0820	<input type="checkbox"/>	ANTHRACITE BLACK	○
CKF0821	<input type="checkbox"/>	SPECIAL DEMAND	€238
SWO4444	<input type="checkbox"/>	DARK RED - SPECIAL	€238
SWO5555	<input type="checkbox"/>	PEARL GREY - SPECIAL	€238
SWO6666	<input type="checkbox"/>	INTENSE YELLOW	€238
SWO7777	<input type="checkbox"/>	VIBRANT ORANGE	€238
SWO8888	<input type="checkbox"/>	POLAR BLUE	€238
SWO9999	<input type="checkbox"/>	CANDY PINK	€238
FRONT CASTORS*			
CKF0850	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID	◆
CKF0870	<input type="checkbox"/>	8" X 2" BLACK SOLID	€32
REAR WHEEL SIZE*			
CKF0942	<input type="checkbox"/>	12" (NC 1111)	◆
REAR WHEEL*			
CKF0963	<input type="checkbox"/>	DRUM BRAKE 4 PARTS FIXING HANDRIM (12" : BLACK PLASTIC HUB AND RIM) DISTANCE WHEEL- HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1120)	○
CKF0964	<input type="checkbox"/>	COMPOSITE 7 PARTS FIXING HANDRIM (BLACK PLASTIC HUB AND RIM) DISTANCE WHEEL- HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1111, 1130)	◆
TYRES*			
CKF0972	<input type="checkbox"/>	BLACK SOLID 1"3/8	◆
REAR WHEEL AXLE*			
CKF1040	<input type="checkbox"/>	FIXED AXLE	◆
REAR WHEEL POSITION*			
CKF1070	<input type="checkbox"/>	PASSIVE RWP	◆
PUMP			
CKF1090	<input type="checkbox"/>	PUMP	€10
CHASSIS OPTION			
CKF1101	<input type="checkbox"/>	FRONT WHL EXTENSION	€30
BRAKES*			
CKF1120	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CKF1130	<input type="checkbox"/>	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)	€243
TRANSPORTATION KIT*			
CKF1175	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆

SAFETY OPTIONS*				
CKF1190		LAP BELT	(NC 1189, 1192, 1193, 1194)	◆
CKF1189		BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON MEDIUM	(NC 1192, 1193, 1194)	€171
CKF1192		BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON LARGE	(NC 1193, 1194)	€171
CKF1193		BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON MEDIUM	(NC 1194)	€161
CKF1194		BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON LARGE		€161
CKF1196		ANTI-TIPPER (PAIR)	(NC 1200, 1245, 1246)	€43
CKF1200		REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(NC 1245, 1246)	€85
OTHER OPTIONS				
CKF1240		CANE HOLDER		€28
CKF1245		TIPPING AID RH		€16
CKF1246		TIPPING AID LH		€16
CKF1270		ACRYLIC TRAY	(NC 1271, 1272, 1273)	€150
CKF1271		HALF TRAY RH		€62
CKF1272		HALF TRAY LH		€62
CKF1273		COMFORT PADDING FOR HALF TRAY		€30
SWO10002		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618		€105

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

CKF0280	NYLON W/O HOOK AND LOOP STRIPS
CKF0295	PADDED NYLON
CKF0300	FIXED 7° BACKREST HT ADJ 15" - 18" FIXED P/HANDLES
CKF0310	FIXED 7° BACKREST HT ADJ 15" - 18" HT ADJ P/HANDLES
CKF0314	FIXED 7° BACKREST 17" ADJ P/HANDLES 6"
CKF0318	ANGLE ADJUSTABLE (-15° TO +15°)
CKF0320	RECLINING (RATCHET) 0/30° HEIGHTS 15"/16"/17"
CKF0334	RECLINING (GAS) 0/30° HEIGHTS 18"/19"/20"
CKF0350	WITHOUT BACK UPHOLSTERY
CKF0360	STANDARD PADDED NYLON SLING TYPE (JEMIMA) - AVAILABLE WITH BACKREST OPTIONS 0300 AND 0310 ONLY
CKF0364	LUMBAR +1 STRAP
CKF0367	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL. BACK POCKET - AVAILABLE WITH BACKREST OPTION 0318 ONLY
CKF0368	CONTOUR 3+1 STRAPS
CKF0375	CONTOUR 7+1 STRAPS
CKF0480	FLIP UP RH
CKF0490	DEPTH AND ANGLE ADJUSTABLE RH
SWO1111	FLIP UP DEPTH AND ANGLE ADJ RH (SEE FULL DESCRIPTION)
CKF0491	ONE PIECE
CKF66666	NO FOOTREST
CKF0500	FLIP UP LH
CKF0510	DEPTH AND ANGLE ADJUSTABLE LH
SWO2222	FLIP UP DEPTH AND ANGLE ADJ LH (SEE FULL DESCRIPTION)
CKF55555	NO FOOTREST
CKF0537	SWING AWAY DETACHABLE RH (N°2) DESK
CKF0538	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH
CKF0539	SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CKF0540	SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CKF0580	REMOVABLE LH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CKF0590	REMOVABLE LH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CKF0593	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CKF0596	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CKF0597	SWING AWAY DETACHABLE LH (N°2) DESK
CKF0598	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH
CKF0599	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CKF0600	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CKF0659	MATRX FLO-TECH LITE CUSHION - Low Risk
CKF0661	MATRX FLO-TECH IMAGE CUSHION - High Risk
CKF0681	MATRX FLO-TECH CONTOUR VISCO NG CUSHION - High Risk
CKF0691	MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk

CKF0701	16" SD
CKF0702	18" SD
CKF0703	20" SD
CKF0656	SEAT RIGIDIZER
CKF0716	NECKREST INCLUDING BACKBRACE
CKF0717	HEADREST INCLUDING BACKBRACE
CKF0718	HEADREST WITH CHEEK SUPPORT
CKF0721	LATERAL POSITIONING PADS PAIR
CKF0760	CALF STRAP
CKF0770	LONGER FOOTREST TUBE
CKF0771	SHORTER FOOTREST TUBE
CKF0780	BACKBRACE
CKF0782	FOLDABLE BACKBRACE
CKF1120	PUSH TO LOCK BRAKES
CKF1130	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)
CKF1189	BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON MEDIUM
CKF1192	BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON LARGE
CKF1193	BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON MEDIUM
CKF1194	BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON LARGE
CKF1200	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CKF1245	TIPPING AID RH
CKF1246	TIPPING AID LH
CKF1270	ACRYLIC TRAY
CKF1271	HALF TRAY RH
CKF1272	HALF TRAY LH
CKF1273	COMFORT PADDING FOR HALF TRAY



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):