



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

BEN NG HD STOCK

Standard Features

- Max user weight 160kg (25 stone)
- Crash tested ISO 7176 - Part 19
- Transit (attendant) propulsion or self (occupant) propulsion
- Seat width range 20" - 22" (510mm - 560mm)
- Seat depth 17" (430mm) or 18" (450mm)
- Seat height 19" (485mm)
- Half folding 10° backrest 16.5" (420mm)
- Black nylon upholstery
- Swing away and detachable legrest hangers with flip up footplates
- Full length, fixed height, removable armrest
- 24" x 1.3/8" quick release rear wheel solid tyre or
- 12" X 2" fixed axle, solid tyre
- 8" x 1.1/4" (200mm x 25mm) solid front castor
- Lapbelt (adjustable both sides)
- Push to lock brakes

Standard Base price: €429



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
Ben NG HD Transit		
1550174	<input type="checkbox"/> 20" W X 17" D	€377
1550175	<input type="checkbox"/> 20" W X 18" D	€377
1550176	<input type="checkbox"/> 22" W X 17" D	€377
1550177	<input type="checkbox"/> 22" W X 18"D	€377
Ben NG Self Propel		
1550184	<input type="checkbox"/> 20" W X 17" D	€429
1550185	<input type="checkbox"/> 20" W X 18" D	€429
1550186	<input type="checkbox"/> 22" W X 17" D	€429
1550187	<input type="checkbox"/> 22" W X 18"D	€429

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
E-mail: ordersireland@invacare.com or Fax: 01 810
7085

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):