



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Ben NG Stock

Standard Features

- Max user weight 127kg (20 stone)
- Crash tested ISO 7176 - Part 19
- Transit (attendant) propulsion or self (occupant) propulsion
- Seat width range 15" - 18" (380mm - 455mm)
- Seat depth 16" (405mm) or 17" (430mm)
- Seat height 19" (485mm)
- Half folding 10° backrest 16.5" (420mm)
- Black nylon upholstery
- Swing away and detachable legrest hangers with flip up footplates
- Full length, fixed height, removable armrest
- 24" x 1.3/8" quick release rear wheel solid tyre or
- 12" X 2" fixed axle, solid tyre
- 8" x 1.1/4" (200mm x 25mm) solid front castor
- Lapbelt (adjustable both sides)
- Push to lock brakes



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
Ben NG Transit		
1550169	<input type="checkbox"/> 15" W X 16" D	€219
1550170	<input type="checkbox"/> 16" W X 16"	€219
1550171	<input type="checkbox"/> 17" W X 17"	€219
1550173	<input type="checkbox"/> 18" W X 17"	€219
Ben NG Self Propel		
1550180	<input type="checkbox"/> 15" W X 16" D	€272
1550181	<input type="checkbox"/> 16" W X 16"	€272
1550182	<input type="checkbox"/> 17" W X 17"	€272
1550183	<input type="checkbox"/> 18" W X 17"	€272

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
E-mail: ordersireland@invacare.com or Fax: 01 810
7085

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):