

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JULY 2018

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Dragon Comfort Seating
Ireland
Max User Weight 140Kg (22 Stone)

Standard Features



- Comfort seat unit with 70mm seat cushion and padded back
- Seat Height 450mm (17.75")
- Seat Widths 440mm - 490mm (17" - 19")
- Seat Depths 410mm - 510mm (16" - 20")
- Detachable T-Armrests, height, width and depth adjustable with standard armpad
- Fixed Headrest
- Manually settable seat Tilt and Electric backrest recline (4°- 20°)
- Standard 80° detachable legrest hangers with adapt footplate
- Shark 2 Remote with Standard Bracket
- Blue Metallic Frame colour
- 230mm Puncture Proof PU-foam sealed, Grey castors
- 318mm x 57mm (12.5" x 2.25") Puncture Proof solid inlay Grey drive wheels
- Rear suspension
- 6 km/h (4mph) 2 Pole Motors and 50 Ah Batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock
- Wheel lock / parking brake
- 8A (200 - 250V) Charger



* Configurations may vary from the image shown

Standard Base Price €4,139

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code		Description		Ireland
AMM0000		Dragon Basic Model		€4,139
AMZ0000		Dragon Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
SEATING TYPE				
AMM1166	X	Comfort Seat Unit	◆	
SEAT CUSHION				
AMM1337	X	Seat Cushion 70mm	◆	
BACK CUSHION				
AMM1101	X	Padded Back Cushion	◆	
UPHOLSTERY TYPE				
AMM0278	X	Grey (Dartex)	◆	
SEAT HEIGHT (Without Seat Cushion)				
AMM0362	X	450mm (17.75") (adjustable from 425mm - 500mm)	◆	
SEAT WIDTH				
AMM0374	X	440mm - 490mm (415mm seat cushion) (17" - 19")	◆	
SEAT DEPTH				
AMM0326		410mm (16")	Not available with Electrical Recliner	○
AMM0327		460mm (18")		○
AMM0330		510mm (20")		○

BACKREST HEIGHT (without Seat Cushion)			
AMM0463	<input checked="" type="checkbox"/>	560mm (22")	◆
ARMREST (Pair)			
AMM0606	<input checked="" type="checkbox"/>	Detachable T-Armrest, adjustable in height 290mm - 360mm	◆
ARMPAD (Pair)			
AMM0647		Standard Armpad (325mm)	◆
AMM0649		Extended Armpad (+50mm)	○
SEAT ACCESSORIES			
HEADREST - Adequate headrest option must be selected			
AMM1427		Headrest with cheek support, sideways adjustable	€31
AMM1354		Headrest Fixed (Clematis)	◆
AMM1495		Equivalent non-Invacare Headrest supporting the requirements of ISO7176-19:2001, provided independently from this order	○
SEAT TILT			
AMM1002		Manually settable Tilt	◆
AMM1014		Electrical Tilt (-4.5° to 25.5°)	€105
BACKREST RECLINER			
AMM0953		Pneumatic Recliner RH (5° to 30°)	Not available with RH control and Tray Bracket ○
AMM0954		Pneumatic Recliner LH (5° to 30°)	Not available with LH Control and Tray Bracket ○
AMM0947		Electrical Recliner (4° - 20°)	◆
FOOTREST HANGER TYPE (Pair)			
AMM0526		Standard Hanger 80°	◆
AMM0529		Manual Height Adjustable (Elevating Legrest)	€42
FOOTPLATE TYPE (Pair)			
AMM0513	<input checked="" type="checkbox"/>	Angle and depth adjustable Footplate (Adapt)	◆
REMOTE CONTROL			
STANDARD REMOTE SHARK 2 & BRACKET (Please also select bracketry option)			
AMM1801	<input checked="" type="checkbox"/>	Shark 2 Remote (enables operation of up to 2 electrical adjustments)	◆
AMM1207		Standard Bracket RH	○
AMM1208		Standard Bracket LH	○
AMM1707		Swing Away + height adjustable RH high end	○
AMM1708		Swing Away + height adjustable LH high end	○
DUAL CONTROL (can only be selected as an addition to user remote)			
AMM1891		Dual Control RH (price applies to dual control and bracket)	€78
AMM1892		Dual Control LH (price applies to dual control and bracket)	€78
DRIVE UNIT OPTIONS			
FRAME COLOUR			
AMM0259	<input checked="" type="checkbox"/>	Blue Metallic (Rear wheel rim Black and front fork Silver)	◆
FRONT WHEEL			
AMM0811		230mm (9") Grey Pneumatic	○
AMM0813		230mm (9") Grey Puncture Proof (Solid Inlay)	◆
REAR WHEEL			
AMM0711		318mm x 57mm (12.5" x 2.25") Grey Pneumatic	○
AMM0727		318mm x 57mm (12.5" x 2.25") Grey Puncture Proof (Solid Inlay)	◆
PERFORMANCE OPTIONS			
SUSPENSION			
AMM1264	<input checked="" type="checkbox"/>	Rear Suspension	◆
MOTOR TYPE			
AMM1282	<input checked="" type="checkbox"/>	6 km/h (4mph) Standard version	◆
BATTERY			
AMM1248		50 Ah (AGM)	◆
AMM1259		60 Ah (Gel)	€42
PERFORMANCE ACCESSORIES			
CHARGER			
AMM1343		OptiCharge Charger 8A (200 - 250V)	◆
AMM1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€150
SAFETY OPTIONS			
ANTI TIPPERS			
AMM1205	<input checked="" type="checkbox"/>	Anti tippers	◆
LIGHTS			
AMM1254	<input checked="" type="checkbox"/>	Without Light system	◆
KERB RIDER			
AMM1244		Central Kerb Rider with front fixation	€210
TRANSPORTATION KIT			
AMM1494	<input checked="" type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory)	◆

FIXATION BELT			
AMM1329		Lap Belt with Auto Lock	◆
ACCESSORIES			
LUGGAGE RACK			
AMM1394	<input checked="" type="checkbox"/>	Rear Carrier	◆
CANE HOLDER			
AMM1312		Cane Holder RH	€42
AMM1363		Cane Holder LH	€42
TRAY - mounted opposite side to controller			
AMM1404		Swing Away RH	€21
AMM1405		Swing Away LH	€21
			Total Price



Special Instructions - Please supply the following information:

Name: _____ **Email:** _____

Phone: _____ **Fax:** _____

Description of special request (please include relevant option code or part number):